

RESPIRATORY CATALOG

EVERYDAY AND SPECIALTY PRODUCTS



Order online at: www.rcmedical.com



FAST DELIVERY

RC MEDICAL, INC.

Phone: 888-431-2727 • Fax: 888-337-0443

2018/2019

P.O. Box 833

• 284 Merrow Road •

Tolland, CT 06084

HEALTHCARE FACILITIES HAVE AN OPEN 30 DAY ACCOUNT

NO MINIMUM ORDER - PURCHASE ORDERS ACCEPTED

F.O.B. Origin. • Prices subject to change without notice.
 25% restocking charge prior to 30 days. • Please call for a RGA Number.
 Products may NOT be returned after 30 days.

INDEX

PAGE	PRODUCT	PAGE	PRODUCT
2-4	Suction Regulators	24	Y-Blocks
5-7	Flowmeters	25	Male Quick Connects
8-9	Respirometers/Repairs	26	Female Couplers
10-12	Cylinder Carts and Stands	27	DISS Fittings
13-16	Regulators	28	O2 Sensors, Analyzers
17	Rail Systems	29	Plastic Bags, Posey Cufflator
18-19	Hoses	30	Compressors, Blenders
20	Hudson Adapters	31	Regulators, Flowmeter Accessories
21-22	Disposables	32	Regulator SPECIALS
23	NPT Fittings		

SUCTION REGULATOR SPECIALS LIMITED QUANTITY



Gentec
Continuous/Intermittent
Suction Regulator
(0-160 mm Hg)



Ohio Medical
"Push to Set" Continuous
2 Mode Suction Regulator



Amvex Digital
Continuous/Intermittent
Suction Regulator
(0 -300 mm Hg)



Amvex Analog
Continuous/Intermittent
Suction Regulator
(0 -300 mm Hg)

	Part #	Price	Part #	Price	Part #	Price	Part #	Price
Gentec	LQ-9	399.00	LQ-10	399.00	LQ-11	399.00	LQ-12	399.00
Ohmeda	LQ-13	279.00	LQ-14	279.00	LQ-15	279.00	LQ-16	279.00
Amvex Digital	LQ-17	429.00	LQ-18	429.00	LQ-19	429.00	LQ-20	429.00
Amvex Analog	LQ-21	399.00	LQ-22	399.00	LQ-23	399.00	LQ-24	399.00



All suction regulators above come standard with tubing nipple outlet.



For DISS Male outlet, add a "D" to end of part number.

Toll Free Phone
888-431-2727
1-860-871-4156

Toll Free Fax
888-337-0443
1-860-870-0443



SUCTION REGULATORS – by Precision Medical CONTINUOUS AND CONTINUOUS/INTERMITTENT



PM3000
2 Mode - Continuous
\$279.00
(no adapters)



PM3100
3 Mode - Continuous
\$279.00
(no adapters)



PM3300
Continuous/Intermittent
\$557.00
(no adapters)

TO ORDER








Use the chart below to determine the inlet and outlet adapters you will need for your facility to obtain the proper part number.

ie: PM3005 - 2 mode regulator with Chemetron Quick Connect and tubing nipple.

Features:


- Compact and lightweight.
- Gauge is fully protected in case.
- Modular design for ease of maintenance.
- Large internal ports assure maximum flows.
- Solid cast metal backplate to ensure strength.
- Large regulator and mode selector knobs for easy adjustment.
- Designed for general patient use with vacuum regulation for 0-200 mm Hg.
- Color coded gauge indicates Low - Medium - High and Full vacuum conditions.
- 10 year warranty



		VACUUM SOURCE CONNECTION				
						
		DISS H.T.	OHMEDA	CHEMETRON	OXEQUIP	PURITAN-BENNETT
PM3000 WITH TUBING NIPPLE		PM3003 308.00	PM3004 314.00	PM3005 314.00	PM3006 314.00	PM3008 311.00
PM3100 WITH TUBING NIPPLE		PM3103 308.00	PM3104 314.00	PM3105 314.00	PM3106 314.00	PM3108 311.00
PM3300 WITH TUBING NIPPLE		PM3303 584.00	PM3304 590.00	PM3305 590.00	PM3306 590.00	PM3308 587.00
PM3000 WITH DISS MALE		PM3014 308.00	PM3015 314.00	PM3016 314.00	PM3017 314.00	PM3019 311.00
PM3100 WITH DISS MALE		PM3114 308.00	PM3115 314.00	PM3116 314.00	PM3117 314.00	PM3119 311.00
PM3300 WITH DISS MALE		PM3314 584.00	PM3315 590.00	PM3316 590.00	PM3317 590.00	PM3319 587.00

VACUTRON SUCTION REGULATORS by Allied

The new, sleekly contoured Vacutron® is a precision engineered instrument that combines the most desirable features for suction regulators. Smaller in size, but big on benefits, it offers advantages that include:




Continuous Suction Regulator

Part #	Inlet	Outlet	Price
22-13-1106	Ohmeda	DISS Male	\$499.00
22-13-1107	DISS Female	DISS Male	499.00
22-13-1108	Chemetron (NCG)	DISS Male	499.00
22-13-1105	Puritan-Bennett	DISS Male	499.00
22-13-1206	Ohmeda	Tubing Nipple	499.00
22-13-1207	DISS Female	Tubing Nipple	499.00
22-13-1208	Chemetron (NCG)	Tubing Nipple	499.00
22-13-1205	Puritan-Bennett	Tubing Nipple	499.00

Continuous/Intermittent Suction Regulator

22-15-1106	Ohmeda	DISS Male	\$899.00
22-15-1107	DISS Female	DISS Male	899.00
22-15-1108	Chemetron (NCG)	DISS Male	899.00
22-15-1105	Puritan-Bennett	DISS Male	899.00
22-15-1206	Ohmeda	Tubing Nipple	899.00
22-15-1207	DISS Female	Tubing Nipple	899.00
22-15-1208	Chemetron (NCG)	Tubing Nipple	899.00
22-15-1205	Puritan-Bennett	Tubing Nipple	899.00

Collection Canisters



Disposable Collection Canisters

Size	Outlet	Units/Case	Part Number	Price
2400ML	DISS	36/CASE	20-08-0001	\$244.00
2400ML	STEM	36/CASE	20-08-0002	237.00
1500ML	STEM	48/CASE	20-08-0003	256.00
1500ML	DISS	48/CASE	20-08-0004	263.00
2400ML	DISS	36/CASE w/tubing	20-08-0005	318.00
2400ML	STEM	36/CASE w/tubing	20-08-0006	337.00
1500ML	STEM	42/CASE w/tubing	20-08-0007	372.00
1500ML	DISS	42/CASE w/tubing	20-08-0008	359.00

Suction Canister Holder and Wall Bracket



White Plastic Holder
20-02-0150
\$19.00



Wall Bracket
22-11-0006
\$49.00

For suction accessories please visit our website at:

www.rcmedical.com

Tubing
Fittings
Vac Trap

Stands
Bacteria Filters
Wye Blocks

PRECISION MEDICAL COMPACT 15 LPM FLOWMETERS



(8MFA1001)

15 LPM OXYGEN FLOWMETER WITH ADAPTER

Part #	Adapter	Price
8MFA1001	None - 1/8" NPT Female	\$36.00
8MFA1002	Diss Female	42.00
8MFA1003	Diss Female w/hand nut	43.00
8MFA1006	Chemetron (NCG)	67.00
8MFA1005	Ohio (Ohmeda)	69.00
8MFA1008	Puritan	55.00
8MFA100M	Medstar	94.00
8MFA1007	Oxequip	52.00

15 LPM OXYGEN FLOWMETER WITH 50 PSI AUXILIARY VALVE AND ADAPTER

Part #	Adapter	Price
8MFA1001 PT0	None - 1/8" NPT Female	\$46.00
8MFA1002 PT0	Diss Female	51.00
8MFA1003 PT0	Diss Female w/hand nut	53.00
8MFA1006 PT0	Chemetron (NCG)	74.00
8MFA1005 PT0	Ohio (Ohmeda)	74.00
8MFA1008 PT0	Puritan	65.00
8MFA100M PT0	Medstar	103.00
8MFA1007 PT0	Oxequip	62.00

15 LPM AIR FLOWMETER WITH ADAPTER

Part #	Adapter	Price
8MFA2001	None - 1/8" NPT Female	\$36.00
8MFA2002	Diss Female	45.00
8MFA2003	Diss Female w/hand nut	50.00
8MFA2006	Chemetron (NCG)	70.00
8MFA2005	Ohio (Ohmeda)	70.00
8MFA2008	Puritan	59.00
8MFA200M	Medstar	98.00
8MFA2007	Oxequip	54.00

15 LPM AIR FLOWMETER WITH 50 PSI AUXILIARY VALVE AND ADAPTER

Part #	Adapter	Price
8MFA2001 PT0	None - 1/8" NPT Female	\$50.00
8MFA2002 PT0	Diss Female	56.00
8MFA2003 PT0	Diss Female w/hand nut	63.00
8MFA2006 PT0	Chemetron (NCG)	78.00
8MFA2005 PT0	Ohio (Ohmeda)	78.00
8MFA2008 PT0	Puritan	73.00
8MFA200M PT0	Medstar	111.00
8MFA2007 PT0	Oxequip	68.00

15 LPM DUAL OXYGEN FLOWMETER WITH ADAPTER

Part #	Adapter	Price
Y8MFA1001	None - 1/8" NPT Female	\$83.00
Y8MFA1002	Diss Female	88.00
Y8MFA1003	Diss Female w/hand nut	90.00
Y8MFA1006	Chemetron (NCG)	113.00
Y8MFA1005	Ohio (Ohmeda)	112.00
Y8MFA1008	Puritan	97.00
Y8MFA100M	Medstar	140.00
Y8MFA1007	Oxequip	98.00

15 LPM DUAL AIR FLOWMETER WITH ADAPTER

Part #	Adapter	Price
Y8MFA2001	None - 1/8" NPT Female	\$88.00
Y8MFA2002	Diss Female	93.00
Y8MFA2003	Diss Female w/hand nut	95.00
Y8MFA2006	Chemetron (NCG)	118.00
Y8MFA2005	Ohio (Ohmeda)	117.00
Y8MFA2008	Puritan	102.00
Y8MFA200M	Medstar	145.00
Y8MFA2007	Oxequip	103.00

15 LPM DUAL OXYGEN FLOWMETER WITH (2) 50 PSI AUXILIARY VALVE AND ADAPTER

Part #	Adapter	Price
Y8MFA1001 PT2	None - 1/8" NPT Female	\$ 93.00
Y8MFA1002 PT2	Diss Female	98.00
Y8MFA1003 PT2	Diss Female w/hand nut	99.00
Y8MFA1006 PT2	Chemetron (NCG)	130.00
Y8MFA1005 PT2	Ohio (Ohmeda)	128.00
Y8MFA1008 PT2	Puritan	106.00
Y8MFA100M PT2	Medstar	149.00
Y8MFA1007 PT2	Oxequip	107.00

15 LPM DUAL AIR FLOWMETER WITH (2) 50 PSI AUXILIARY VALVE AND ADAPTER

Part #	Adapter	Price
Y8MFA2001 PT2	None - 1/8" NPT Female	\$ 98.00
Y8MFA2002 PT2	Diss Female	103.00
Y8MFA2003 PT2	Diss Female w/hand nut	104.00
Y8MFA2006 PT2	Chemetron (NCG)	135.00
Y8MFA2005 PT2	Ohio (Ohmeda)	133.00
Y8MFA2008 PT2	Puritan	111.00
Y8MFA200M PT2	Medstar	154.00
Y8MFA2007 PT2	Oxequip	112.00



(Y8MFA1003PT2)

CHAIN-RING ASSEMBLY



The Chain-Ring Assembly is used to permanently attach Flowmeters to the wall outlet.

Part #	Description	Price
8004	6" Chain w/Closed Ring	\$7.00
8005	6" Chain w/Open Ring	\$7.00

ADAPTERS



Diss Female



Diss Female with hand nut



Chemetron (NCG)



Ohio (Ohmeda)



Puritan



Medstar



Oxequip



50 PSI Auxiliary Valve

(oxygen) OPTO - \$16.00

(air) APTO - \$19.00

PRECISION MEDICAL CHROME BODY FLOWMETER

**5 Year
Warranty**



4MFA1001



1MFA9001



3MFA1001



0 - 200 cc PEDIATRIC OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
6MFA1001	None - 1/8" NPT Female	\$269.00
6MFA1003	DISS FEMALE w/Hand Tight	275.00
6MFA1006	CHEMETRON (NCG)	284.00
6MFA1005	OHIO/OHMEDA	285.00
6MFA1008	PURITAN BENNETT	279.00

0 - 1 LPM PEDIATRIC OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
4MFA1001	None - 1/8" NPT Female	\$230.00
4MFA1003	DISS FEMALE w/Hand Tight	243.00
4MFA1006	CHEMETRON (NCG)	248.00
4MFA1005	OHIO/OHMEDA	252.00
4MFA1008	PURITAN BENNETT	243.00

0 - 3 LPM PEDIATRIC OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA3001	None - 1/8" NPT Female	\$96.00
1MFA3003	DISS FEMALE w/Hand Tight	103.00
1MFA3006	CHEMETRON (NCG)	127.00
1MFA3005	OHIO/OHMEDA	127.00
1MFA3008	PURITAN BENNETT	105.00

0 - 8 LPM OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
8MFA4001	None - 1/8" NPT Female	\$36.00
8MFA4003	DISS FEMALE w/Hand Tight	43.00
8MFA4006	CHEMETRON (NCG)	67.00
8MFA4005	OHIO/OHMEDA	67.00
8MFA4008	PURITAN BENNETT	55.00

0 - 70 LPM HI FLOW OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA8001	None - 1/8" NPT Female	\$106.00
1MFA8003	DISS FEMALE w/Hand Tight	112.00
1MFA8006	CHEMETRON (NCG)	124.00
1MFA8005	OHIO/OHMEDA	125.00
1MFA8008	PURITAN BENNETT	117.00

0 - 70 LPM HI FLOW AIR FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA9001	None - 1/8" NPT Female	\$106.00
1MFA9003	DISS FEMALE w/Hand Tight	119.00
1MFA9006	CHEMETRON (NCG)	124.00
1MFA9005	OHIO/OHMEDA	126.00
1MFA9008	PURITAN BENNETT	120.00

SELECT FLOWMETER

15 LPM OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
3MFA1001	None - 1/8" NPT Female	\$80.00
3MFA1003	DISS FEMALE w/Hand Tight	88.00
3MFA1006	CHEMETRON (NCG)	98.00
3MFA1005	OHIO/OHMEDA	100.00
3MFA1008	PURITAN BENNETT	92.00

FOR MORE FLOWMETERS, PLEASE VISIT OUR WEBSITE AT:
WWW.RCMEDICAL.COM


FLOW SELECTOR



Flow Selector Features

- Give nebulizer treatments without disconnecting the humidifier
- Reduces the changeover time to seconds
- No more loose tubing laying around
- Eliminates tubing nipple adaptors (It has 2 of its own.)
- Fits all flowmeters that have DISS male outlet fittings
- Lightweight

Part # 1000 Price \$41.00




INCREMENT SCALES

0-200cc, 0-1LPM, 0-3 LPM,
0-70 LPM, 0-15 LPM

FIXED FLOW OXYGEN FLOWMETERS (1/2 - 15 LPM)

Position Increments (LPM)
Off, .50, 1, 2, 3, 4, 5, 6, 8, 10, 15



It's Durable - The lack of a flow tube and the durable metal body prevents breakage problems long associated with standard flowmeters.

It's Compact - The compact size allows the unit to fit into tight headwall areas and makes it virtually damage proof when other equipment is installed.





It's Accurate - The patented multi-spring technology insures accuracy at every flow setting without the need for a visible float as in Thorpe tube flowmeters

It's User-Friendly - The knurled flow adjustment knob and the large, easy-to-read numbers make flow adjustment easy.

It's Versatile - The unit is available for headwall outlets as well as for ceiling mount drops. Both air and oxygen models are available.

FOR MORE FLOWMETERS, PLEASE VISIT
OUR WEBSITE AT:
WWW.RCMEDICAL.COM



	Adapters	Side Connection	Top Connection
	1/8" NPT FE	7MFA1001 \$65.00	7MFA1001TC \$65.00
	DISS Female Hand Tight	7MFA1003 \$71.00	7MFA1003TC \$71.00
	Ohmeda	7MFA1005 \$90.00	7MFA1005TC \$90.00
	Chemtron (NCG)	7MFA1006 \$88.00	7MFA1006TC \$88.00
	Puritan-Bennett	7MFA1008 \$84.00	7MFA1008TC \$84.00




NEB BLOCK w/FLOWMETER

Features:

- The Neb Block allows a gas outlet to simultaneously or individually provide various respiration treatments without disconnecting any tubing or equipment.
- Provides easy to use On/Off switch to quickly access oxygen flow to a hand held nebulizer or resuscitation bag.
- Reliable fixed flow orifices provide six ranges of flow for aerosol therapy or resuscitation.



6-7 LPM 8-9 LPM 14-16 LPM
7-8 LPM 10-11 LPM

	6-7 LPM	8MFA6065	137.00
	7-8 LPM	8MFA6075	137.00
	8-9 LPM	8MFA6085	137.00
	10-11 LPM	8MFA6105	137.00
	14-16 LPM	8MFA6155	137.00
	6-7 LPM	8MFA6066	136.00
	7-8 LPM	8MFA6076	136.00
	8-9 LPM	8MFA6086	136.00
	10-11 LPM	8MFA6106	136.00
	14-16 LPM	8MFA6156	136.00
	6-7 LPM	8MFA6063	136.00
	7-8 LPM	8MFA6073	136.00
	8-9 LPM	8MFA6083	136.00
	10-11 LPM	8MFA6103	136.00
	14-16 LPM	8MFA6153	136.00
	6-7 LPM	8MFA6068	130.00
	7-8 LPM	8MFA6078	130.00
	8-9 LPM	8MFA6088	130.00
	10-11 LPM	8MFA6108	130.00
	14-16 LPM	8MFA6158	130.00

SPECIALTY OXYGEN FLOWMETERS



Precision Medical 15 LPM O2 Flowmeter for Ohmeda Wall Outlet with 50 PSI-Ohmeda Female Quick Connect
SFO-1 \$189.00



Precision Medical 15 LPM O2 Flowmeter for Chemetron Wall Outlet with 50 PSI-Chemetron Female Quick Connect
SFO-2 \$169.00



Precision Medical 15 LPM O2 Flowmeter for Puritan-Bennett Wall Outlet with 50 PSI-Puritan-Bennett Female Quick Connect
SFO-3 \$169.00



Precision Medical 15 LPM O2 Flowmeter for Chemetron Wall Outlet with 50 PSI-Chemetron Female Quick Connect
SFO-4 \$217.00



Precision Medical 15 LPM O2 Flowmeter for Ohmeda Wall Outlet with 50 PSI-Ohmeda Female Quick Connect
SFO-5 \$217.00



Precision Medical 15 LPM O2 Flowmeter for Puritan-Bennett Wall Outlet with 50 PSI-Puritan Bennett Female Quick Connect
SFO-6 \$217.00

RESPIROMETERS

MARK 8



1.4" Face

For measurement of minute volumes as well as tidal volumes. Two dials. Small dial records 0-1 liter. Large dial 0-100 liters.

MARK 14



1.4" Face

Specially suited for the measurement of tidal volumes and may also be used for measurement of minute volumes. Two dials. Small dial records 0-10 liters. Large dial records 0-1 liter.

HALOSCALE STANDARD



2.25" Face

The original design of Haloscale Respirometer incorporating concentric dials and an enlarged face for improved readability. Equally convenient for measuring both tidal and minute volumes. Large 57 mm dial diameter. Red/black hands and black scales on a white background with yellow highlight ring. Two concentric outer dials. First dial records 0-1 liter. Second dial records 0-100 liters.

ORDER #	RESPIROMETERS DESCRIPTION	PRICE
2750.....	Mark 8 Respirometer	\$2105.00
2751.....	Mark 14 Respirometer	2105.00
2768.....	Haloscale Standard.....	2155.00

RESPIROMETERS WITH SPIRO GUARD

2756.....	Mark 8 Respirometer	2169.00
2757.....	Mark 14 Respirometer	2169.00
2769.....	Haloscale Standard.....	2219.00

RESPIROMETERS WITH SPIRO GUARD AND EXTRA LARGE CARRYING CASE

2762.....	Mark 8 Respirometer	2199.00
2763.....	Mark 14 Respirometer	2199.00
2772.....	Haloscale Standard.....	2249.00

Used Respirometers

Fully Serviced, Cleaned,
Oiled and Calibrated

1 Year Warranty

Order #2799 \$1150
(Case included)

ELECTRONIC DIGITAL RESPIROMETER

A respirometer long overdue in this modern electronic age




ANESTHESIA ASSOCIATES

MODEL 295 DIGITAL RESPIROMETER IS AN ELECTRONIC VERSION OF THE MECHANICAL RESPIROMETER. IT DOES THE SAME JOB PLUS

- Easier to read
- Saving on the initial purchase and maintenance
- More reliable
- Extremely durable
- No gears, shafts, or bearings to break down in the electronic unit
- Separates for ease of patient use, sterilization and cleaning

Part #		Price
295	Respirometer	\$1190.00
Options		
296	6' Retractable Cable (extended length)	\$144.00
298	Carrying Case	88.00
299	Spare Flow Sensing Unit	395.00

Model 295, a two piece detachable respirometer comes complete ready to use with 9v battery and one (1) 22 mm female x 22 mm female adapter.



All respirometers come with a standard zippered carrying case, two metal adapters, and harness.

Spiro-Guards and extra large carrying cases are optional.

RUBBER RESPIROMETER PROTECTORS

2700 Spiro-Guard



2701 Spiro-Guard



2702 Spiro-Guard



Minimizes damage caused by impact
Protects vulnerable external parts
One piece 1/4 inch natural rubber covers entire outer body
Allows full view of dial as well as access to On/Off and Reset buttons

2700 Designed to fit: mark 8, mark 14, mark 25, physiological, haloscale compact, wright, penlon\$86.00
 2701 Designed to fit: All Fraser-Harlake models OHMEDA (3 buttons).....96.00
 2702 Designed to fit: mark 12, mark 20, haloscale standard, wright86.00

CARRYING CASES

REPAIRS


2703 Carrying Case\$67.00



The 2703 carrying case was specifically designed EXTRA LARGE for the purpose of holding respirometers with Spiro-Guards on them. The plastic case does NOT have a zipper and is lightweight and durable.

(Respirometers, Spiro-Guards, and adapters not included.) 9" x 5" with carrying handle

Replacement Standard Zippered Carrying Case
 #108517.....\$79.00



Respirometers - \$199.00
 includes cleaning, calibration and oiling. Parts are extra.

It is not necessary to obtain an RGA to send in repairs.

Please send them to the address below. Most repairs are returned to you in 5 days.

RC Medical Inc.
 284 Merrow Road
 Tolland, CT 06084

In our 20 years of respirometer repairing, these are the most common problems, in addition to our recommended yearly cleaning:

- Lens Cover cracked or missing
- Vane bent or blown out
- Rust

Most repairs are \$199 - \$275

All Models repaired including:

Wright • Ferraris
Haloscale • BOC
Fraser - Harlake
Ohmeda
Boehringer

WE WILL GIVE YOU UP TO \$400.00
IN FREE GOODS, FOR EVERY RESPIROMETER YOU
SEND IN FOR TRADE.



SMALL CYLINDER CARTS (D&E)

2-Wheel Cart



Part #	Cylinders Held	Price	Dimensions (Height-Depth-Width)
6041	4	181.00	42" - 10" - 18"
6061	6	189.00	42" - 15" - 18"
6121	12	203.00	42" - 20" - 22"

With Sound Deadening Kit (below)

6041-SDK	4	214.00	42" - 10" - 18"
6061-SDK	6	233.00	42" - 15" - 18"
6121-SDK	12	247.00	42" - 20" - 22"

4-Wheel Cart



Part #	Cylinders Held	Price	Dimensions (Height-Depth-Width)
6064	6	254.00	42" - 15" - 18"
6124	12	274.00	42" - 20" - 22"
6184*	18	282.00	42" - 30" - 16"
6244*	24	303.00	42" - 30" - 20"

With Sound Deadening Kit (below)

6124-SDK	12	311.00	42" - 20" - 22"
6184-SDK*	18	323.00	42" - 30" - 16"
6244-SDK*	24	379.00	42" - 30" - 20"

* -4" wheels

Multiple D & E Size Cylinder Stands With Wheels



Part #	Cylinders Held	Price	Dimensions (Height-Depth-Width)
6040-PC	4	164.00	23" - 13" - 10"
6060-PC	6	181.00	23" - 12" - 15"
6120-PC	12	198.00	23" - 15" - 22"

With Sound Deadening Kit (below)

6040-PC-SDK	4	193.00	23" - 13" - 11"
6060-PC-SDK	6	210.00	23" - 12" - 15"
6120-PC-SDK	12	223.00	23" - 15" - 22"

ALL CARTS ABOVE ARE PAINTED GREEN

RED CARTS - ADD A "R" TO END OF PART NUMBER, PRICE IS \$40.00 MORE
YELLOW CARTS - ADD A "Y" TO END OF PART NUMBER, PRICE IS \$40.00 MORE

SOUND DEADENING KITS (for Carts Above)



Anthony has developed an exclusive sound deadening system for their multiple cylinder carts, used to distribute cylinders in and around hospitals, rest homes, convalescent homes, etc. or where QUIET is essential.

These sound deadening kits include heavy-duty sound deadening matting and plastic bumper guards. This formula eliminates the rattling & clanging of steel & aluminum cylinders while in transport.

SOUND DEADENING KITS ALSO SOLD SEPARATELY
PLEASE VISIT OUR WEBSITE: WWW.RCMEDICAL.COM



10" x 14" Plastic Sign
Part # SP-51
\$37.00



10" x 14" Plastic Sign
Part # SP-52
\$37.00

SMALL CYLINDER STANDS AND HOLDERS (D&E)

Multiple D & E Size Cylinder Stands



Part #	Cylinders Held	Price	Dimensions (Height-Depth-Width)
6040	4	143.00	20" - 13" - 11"
6060	6	153.00	20" - 12" - 15"
6120	12	163.00	20" - 15" - 22"
6180	18	180.00	20" - 15" - 32"
6240	24	195.00	20" - 20" - 32"
6400	40	315.00	20" - 40" - 25"

With Sound Deadening Kit (below)

6040-SDK	4	172.00	20" - 13" - 11"
6060-SDK	6	187.00	20" - 12" - 15"
6120-SDK	12	201.00	20" - 15" - 22"
6180-SDK	18	232.00	20" - 15" - 32"
6240-SDK	24	264.00	20" - 20" - 32"
6400-SDK	40	421.00	20" - 40" - 25"

Can be bolted to floor or free standing.

ALL STANDS ABOVE ARE PAINTED GREEN

RED STANDS - ADD A "R" TO END OF PART NUMBER, PRICE IS \$40.00 MORE

YELLOW STANDS - ADD A "Y" TO END OF PART NUMBER, PRICE IS \$40.00 MORE



MULTIPLE CYLINDER STANDS

CHROME



PX-1025

PX-1025 "D/E" Rack / 6 cylinders 9" x 15" x 10" \$ 93.00

PX-1027 "D/E" Rack / 12 cylinders 9" x 15" x 20" 145.00

Bolts to floor or free standing.



For any solid surfaces,
crash carts, more...

Part # 610WM
\$67.00



For wheelchairs
Part # 610WC
\$77.00

FOR MORE INFORMATION ON SMALL CYLINDER BASES,
CYLINDER HOLDERS FOR BEDS, IV POLES, M6 CYLINDER CARTS AND STANDS,
AND MORE, PLEASE VISIT OUR WEBSITE AT:

WWW.RCMEDICAL.COM

SMALL CYLINDER CARTS (D&E)

#6105.....\$49.00



#PX-1060.....\$49.00



#6103-4.....\$173.00



These four-wheeled units offer better control and easy maneuverability in and out of tight spaces. They require no tilting or balancing, giving you total control over the cylinder load.

#6106-CH
\$79.00

Cane Handle



#13001-2
\$112.00

Dual Tank



LARGE CYLINDER CARTS AND STANDS (H,K)

Dual
6214



Single
6114



All models have 4 wheels.

Part #	Description	Price
6114	Single (Most Popular)	\$287.00
6114-PN	Single (Pneumatic Wheels)	327.00
6214	Dual	339.00
6214-PN	Dual (Pneumatic Wheels)	379.00

2 Wheels Only



Part # HCC2-M
\$187.00




610FW-E.....\$189



Fiberglass Cylinder Bases
HS-2.....\$73

HELIOX, AIR, CO₂, NITROGEN REGULATORS

Heliox Regulators



0 - 20 LPM

Part #	Ratio	For Cylinder size	Price
2210-20	70 - 30	Small D and E CGA 890 Yoke	\$297
2211-20	80 - 20	Small D and E CGA 890 Yoke	\$297
2212-20	70 - 30	Large CGA 280 Nut/Nipple	\$289
2213-20	80 - 20	Large CGA 280 Nut/Nipple	\$289

Nitrogen Regulator with Schrader Coupler















CGA 580 large Cylinder

PSI	Part #	Price
0-180	SRO-45	\$489.00
0-450	SRO-46	\$489.00

WESTERN MEDICA MEDICAL REGULATORS (AIR or CARBON DIOXIDE)


Single-stage regulators preset to 50 psi or adjustable to 100 psi, or with a flow-meter. A sensitive neoprene diaphragm assembly provides optimum accuracy and sensitivity. Dual inlet filtering. Available for small or large cylinders.

	AIR	CO ₂
Preset 50 PSI	 M1A-346P \$188.00 for large cylinder	 M1-320P \$180.00 for large cylinder
Adjustable 0 - 100 PSI	 M1-346PG \$204.00 for large cylinder	 M1-320PG \$208.00 for large cylinder
Liters Per Minute (lpm)	 M1-346-15FM \$232.00 for large cylinder	 M1-320-12FM \$230.00 for large cylinder
Preset 50 PSI	 M1A-950P \$193.00 for small cylinder	 M1-940P \$193.00 for small cylinder
Adjustable 0 - 100 PSI	 M1-950PG \$209.00 for small cylinder	 M1-940PG \$208.00 for small cylinder
Liters Per Minute (lpm)	 M1-950-15FM \$229.00 for small cylinder	 M1-940-12FM \$258.00 for small cylinder

DUAL NITROGEN REGULATOR

Dual nitrogen single-stage high purity regulator. Delivery Range: 0 - 180 psi. 6 year warranty. Comes with Schrader quick connects as shown.

Part # SRO-47.....\$887.00



FOR MORE REGULATORS
PLEASE VISIT OUR WEBSITE
WWW.RCMEDICAL.COM

OXYGEN REGULATORS - LITERS PER MINUTE

Allied All Brass Oxygen Regulators

0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 settings



Part #21017.....\$109.00
for small cylinder (870)

Responsive Respiratory Fixed Flow E Cylinder Oxygen Regulator

0, .12, .5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15 settings




Part #120-3040\$35.00
for small cylinder (870)



Part #21025.....\$169.00
1 - 25 LPM

Two High Pressure Outlets Also

0, .12, .5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15 settings



Part #120-3060\$35.00
for large cylinder (540)

Western Medica Oxygen Regulators

(Diaphragm Type)

6 Year Warranty



M1-540-15 FM

Single stage flowmeter/regulator for large cylinders.
15 liters per minute.
ORDER #M1-540-15 FM.....\$191.00

Single stage flowmeter/regulator for small cylinders.
15 liters per minute.
ORDER #M1-870-15 FM.....\$197.00



15 LPM

Single stage, swivel yoke regulator for small cylinders.
15 liters per minute.
ORDER #MR-870-15 FG\$149.00

With Gauge-Savers
ORDER #MR-870-15 FGP\$173.00

All Western Medica Oxygen Regulators include plastic nut & nipple



15 LPM

Single stage regulator for large cylinders.
15 liters per minute.
ORDER #M1-540-15 FG\$149.00

With Gauge-Savers
ORDER #M1-540-15 FGP\$173.00



15 LPM


Single stage regulator for small cylinders.
15 liters per minute.
ORDER #M1-870-15 FG\$149.00

With Gauge-Savers
ORDER #M1-870-15 FGP\$173.00

To order 8 or 5 liter regulators, add an 8 or a 5 to the end of our part number. The same price as our 15 liter regulators.



Gauge Savers (Rubber)



MADE IN USA



Gauge Savers Help:

- reduce repair costs
- save downtime
- keep regulators "in service"

FOR MORE
OXYGEN REGULATORS
PLEASE VISIT OUR
WEBSITE AT:
WWW.RCMEDICAL.COM

HIGH PRESSURE REGULATORS



Adjustable 0 - 100 psi with DISS Outlet

Gas	Part #	Cylinder Size	Price
Oxy.	MI 540 PG	CGA 540 (Large)	\$167.00
Oxy.	MI 870 PG	CGA 870 (Small)	\$154.00
Air	MI 346 PG	CGA 346 (Large)	\$204.00
Air	MI 950 PG	CGA 950 (Small)	\$209.00



Preset 50 psi with DISS Outlet

Gas	Part #	Cylinder Size	Price
Oxy.	MIA 540 P	CGA 540 (Large)	\$154.00
Air	MIA 346 P	CGA 346 (Large)	\$188.00
N ₂ O	MIA 326 P	CGA 326 (Large)	\$154.00
CO ₂	MIA 320 P	CGA 320 (Large)	\$180.00



Preset 50 psi with DISS Outlet

Gas	Part #	Cylinder Size	Price
Oxy.	MIA 870 P	CGA 870 (Small D&E)	\$131.00
Air	MIA 950 P	CGA 950 (Small D&E)	\$193.00
N ₂ O	MIA 910 P	CGA 910 (Small D&E)	\$158.00
CO ₂	MIA 940 P	CGA 940 (Small D&E)	\$193.00

MADE IN USA



Preset 50 psi with Ohmeda Quick Connect

Gas	Part #	Cylinder Size	Price
Oxy.	SRO-33	CGA 540 (Large)	\$232.00
Air	SRO-34	CGA 346 (Large)	\$252.00



Preset 50 psi with Chemetron Quick Connect

Gas	Part #	Cylinder Size	Price
Oxy.	SRO-35	CGA 540 (Large)	\$232.00
Air	SRO-36	CGA 346 (Large)	\$252.00



Preset 50 psi with Ohmeda Quick Connect

Gas	Part #	Cylinder Size	Price
Oxy.	SRO-30	CGA 870 (Small D&E)	\$234.00
Air	SRO-32	CGA 950 (Small D&E)	\$254.00



Preset 50 psi with Chemetron Quick Connect

Gas	Part #	Cylinder Size	Price
Oxy.	SRO-29	CGA 870 (Small D&E)	\$234.00
Air	SRO-31	CGA 950 (Small D&E)	\$254.00

SPECIALTY OXYGEN REGULATORS



Allied Healthcare Fixed Flow Regulator with Ohmeda 50 PSI Outlet and DISS 50 PSI Outlet, 0 - 25 LPM
Part # SRO-11..... \$274.00



Allied Healthcare Fixed Flow Regulator with Chemetron 50 PSI Outlet and DISS 50 PSI Outlet, 0 - 25 LPM
Part # SRO-12..... \$274.00



Allied Healthcare Fixed Flow Regulator with Puritan Bennett 50 PSI Outlet and DISS 50 PSI Outlet, 0 - 25 LPM
Part # SRO-13..... \$274.00



Oxygen Regulator/Flowmeter 15 LPM, Large Cylinder with 50 PSI Valve
Part # SRO-14..... \$274.00



Oxygen Regulator with Dual Chemetron Female Couplers, 50 PSI, Large Cylinder
Part # SRO-15..... \$304.00



Oxygen Regulator with Dual Ohmeda Female Couplers, 50 PSI, Large Cylinder
Part # SRO-16..... \$304.00



Dual Oxygen Regulator, Large Cylinder
Part # SRO-17..... \$279.00



Dual Oxygen Regulator with (2) 50 PSI Auxiliary Valves, Large Cylinder
Part # SRO-18..... \$319.00



Oxygen Regulator, Preset 50 PSI with Ohmeda Quick Connect, Large Cylinder
Part # SRO-19..... \$232.00



Oxygen Regulator, Preset 50 PSI with Puritan Bennett Quick Connect, Large Cylinder
Part # SRO-20..... \$222.00



Oxygen Regulator, Preset 50 PSI with Chemetron Quick Connect, Large Cylinder
Part # SRO-21..... \$232.00



Oxygen Regulator with 15 LPM flowmeter and 50 PSI Ohmeda Quick Connect, Large Cylinder
Part # SRO-22..... \$329.00



Oxygen Regulator with 15 LPM Flowmeter and 50 PSI Chemetron Quick Connect, Large Cylinder
Part # SRO-23..... \$329.00



Flotec 0 - 25 LPM Regulator with Puritan Bennett Female 50 PSI Outlet for Large Cylinder (CGA 540)
Part # SRO-24..... \$159.00



Flotec 0 - 25 LPM Regulator with Chemetron Female 50 PSI Outlet for Large Cylinder (CGA 540)
Part # SRO-25..... \$159.00




Flotec 0 - 25 LPM Regulator with Ohmeda Female 50 PSI Outlet for Large Cylinder (CGA 540)
Part # SRO-26..... \$159.00

Specialty AIR Regulators, Visit our Website at: WWW.RCMEDICAL.COM

For Regulator Accessories See Page 31

EASTERN RAIL SYSTEMS

Simple Outlet Expansion allows your hospital to use its existing outlets, and to create the number and placement of outlets needed. Simply mount a length of our rail on the wall - attach the needed outlets to the rail with our thumb screw clamp - attach a length of high-pressure hose to the block, and run it back to your existing outlet.



- Easy external mounting - no need to shut down the unit.
- No more weight on outlets
- Shift the weight of humidifiers & suction bottles to the rail
- Eliminate overloading
- Reduce outlet leakage. Save time and money
- Increase efficiency
- Works with all piping systems
- Have outlets where you need them NOW!


Kit comes with:

- 1 each 2 Foot Rail
- 2 each Mounting Hardware
- 1 each 3 Foot Hose with Wall Adapter
- 1 each Double Gas Block

Outlet connection	GAS	Part #	Price
DISS	Oxygen	ERO-1DS	479.00
Chemetron	Oxygen	ERO-1CH	597.00
Ohmeda	Oxygen	ERO-10M	597.00
DISS	AIR	ERA-1DS	479.00
Chemetron	AIR	ERA-1CH	597.00
Ohmeda	AIR	ERA-10M	597.00
DISS	VAC	ERV-1DS	479.00
Chemetron	VAC	ERV-1CH	597.00
Ohmeda	VAC	ERV-10M	597.00



**FOR MORE EASTERN RAIL KITS
AND ACCESSORIES,
PLEASE GO TO OUR WEBSITE AT:
WWW.RCMEDICAL.COM**



CUSTOM HOSE ASSEMBLIES

Ordering is easy as 1, 2, 3 . . . 4


How to Order:

Western Medica manufactures custom hose assemblies that conform to CGS, and NFPA standards. All assemblies are pressure-tested. Select your requirement from the charts below.

Example: OXY-024-FNN-OHF

- 1) OXY Oxygen Service
- 2) 024 24" Long
- 3) FNN Inlet Connection
- 4) OHF Outlet Connection

Evacuation Hose (purple) price may be more, depending on your adaptor



(Pricing for the above example)

FNN	17.00
OHF	69.00
Hose (.35 per inch x 24")	8.40
<u>Assembly Charge</u>	<u>12.00</u>
	\$106.40

MINIMUM HOSE LENGTH — 1 FOOT

Gas Service

Gas	Service	Color
OXY	Oxygen	Green
AIR	Medical Air	Yellow
SUC	Suction	White
N2O	Nitrous Oxide	Blue
EVC	Evacuation	Purple
CO ₂	Carbon Dioxide	Gray
NIT	Nitrogen	Black
GNR	General	Clear

Adapters pictured below may be ordered Separately.

When ordering adapters only, please specify gas.

Example: chm - oxy

All hose assemblies are 1/4" conductive (unless specified)

 #CHM \$38.00 Chemetron (NCG) Male	 #FNN \$17.00 Diss Female with Hex nut	 #FHT \$19.00 Diss Female with Hand Nut	 #OHM \$39.00 Ohmeda Male
 #MBA \$19.00 Diss Male	 #SHM \$21.00 Schrader Male	 #OXM \$24.00 Oxequip Male	 #PBM \$24.00 Puritan - Bennett Male
 #CHF \$67.00 Chemetron (NCG) Female	 #SHF \$87.00 Schrader Female	 #PBF \$61.00 Puritan - Bennett Female	 #OHF \$69.00 Ohmeda Female
 #DBA \$20.00 Demand Body Adapter	 #MP2 \$5.00 Male Pipe 1/8"	 #MP4 \$5.00 Male Pipe 1/4"	
 #FP2 \$5.00 Female Pipe 1/8"	 #FP4 \$5.00 Female Pipe 1/4"	 #BLK \$0.00 Blank (No Connection)	

REINFORCED MEDICAL HOSE



Hose only!
\$4.20 per foot

Conductive, Non-UL Listed
400 psi Working Pressure
1500 psi Burst Pressure

MH4-04C.....	1/4" ID, Nitrous Oxide (Blue)
MH4-08C.....	1/4" ID, Carbon Dioxide (Grey)
MH4-12C.....	1/4" ID, Nitrogen (Black)
MH4-16C.....	1/4" ID, Medical Air (Yellow)
MH4-22C.....	1/4" ID, Suction (White)
MH4-24C.....	1/4" ID, Oxygen (Green)
MH5-22C.....	5/16" ID, Suction (White)

FERRULE



For 1/4" hose
.525 ID
#M-625-A.....\$1.00

VICE-GRIPS

For 3/16", 1/4" and
5/16" ID Hose
C-5A \$59.00



COMMON HOSE ASSEMBLIES

DISS FEMALE HANDTIGHT X DISS FEMALE HANDTIGHT



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-100	\$69.00
5 Foot Air	CH-101	\$69.00
10 Foot Oxygen	CH-102	\$90.00
10 Foot Air	CH-103	\$90.00
15 Foot Oxygen	CH-104	\$111.00
15 Foot Air	CH-105	\$111.00
20 Foot Oxygen	CH-106	\$132.00
20 Foot Air	CH-107	\$132.00

DISS FEMALE HANDTIGHT X OHMEDA MALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-132	\$87.00
5 Foot Air	CH-133	\$87.00
10 Foot Oxygen	CH-134	\$108.00
10 Foot Air	CH-135	\$108.00
15 Foot Oxygen	CH-136	\$129.00
15 Foot Air	CH-137	\$129.00
20 Foot Oxygen	CH-138	\$150.00
20 Foot Air	CH-139	\$150.00

DISS FEMALE HANDTIGHT X CHEMETRON MALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-108	\$86.00
5 Foot Air	CH-109	\$86.00
10 Foot Oxygen	CH-110	\$107.00
10 Foot Air	CH-111	\$107.00
15 Foot Oxygen	CH-112	\$128.00
15 Foot Air	CH-113	\$128.00
20 Foot Oxygen	CH-114	\$149.00
20 Foot Air	CH-115	\$149.00

DISS FEMALE HANDTIGHT X OHMEDA FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-140	\$119.00
5 Foot Air	CH-141	\$119.00
10 Foot Oxygen	CH-142	\$140.00
10 Foot Air	CH-143	\$140.00
15 Foot Oxygen	CH-144	\$161.00
15 Foot Air	CH-145	\$161.00
20 Foot Oxygen	CH-146	\$182.00
20 Foot Air	CH-147	\$182.00

DISS FEMALE HANDTIGHT X CHEMETRON FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-116	\$117.00
5 Foot Air	CH-117	\$117.00
10 Foot Oxygen	CH-118	\$138.00
10 Foot Air	CH-119	\$138.00
15 Foot Oxygen	CH-120	\$159.00
15 Foot Air	CH-121	\$159.00
20 Foot Oxygen	CH-122	\$180.00
20 Foot Air	CH-123	\$180.00

OHMEDA MALE X OHMEDA FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-148	\$138.00
5 Foot Air	CH-149	\$138.00
10 Foot Oxygen	CH-150	\$159.00
10 Foot Air	CH-151	\$159.00
15 Foot Oxygen	CH-152	\$180.00
15 Foot Air	CH-153	\$180.00
20 Foot Oxygen	CH-154	\$201.00
20 Foot Air	CH-155	\$201.00

CHEMETRON MALE X CHEMETRON FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-124	\$135.00
5 Foot Air	CH-125	\$135.00
10 Foot Oxygen	CH-126	\$156.00
10 Foot Air	CH-127	\$156.00
15 Foot Oxygen	CH-128	\$177.00
15 Foot Air	CH-129	\$177.00
20 Foot Oxygen	CH-130	\$198.00
20 Foot Air	CH-131	\$198.00

PURITAN BENNETT MALE X PURITAN BENNETT FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-156	\$115.00
5 Foot Air	CH-157	\$115.00
10 Foot Oxygen	CH-158	\$136.00
10 Foot Air	CH-159	\$136.00
15 Foot Oxygen	CH-160	\$157.00
15 Foot Air	CH-161	\$157.00
20 Foot Oxygen	CH-162	\$178.00
20 Foot Air	CH-163	\$178.00

DISS FEMALE HANDTIGHT X PURITAN BENNETT MALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-164	\$72.00
5 Foot Air	CH-165	\$72.00
10 Foot Oxygen	CH-166	\$93.00
10 Foot Air	CH-167	\$93.00
15 Foot Oxygen	CH-168	\$114.00
15 Foot Air	CH-169	\$114.00
20 Foot Oxygen	CH-170	\$135.00
20 Foot Air	CH-171	\$135.00

DISS FEMALE HANDTIGHT X PURITAN BENNETT FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-172	\$111.00
5 Foot Air	CH-173	\$111.00
10 Foot Oxygen	CH-174	\$132.00
10 Foot Air	CH-175	\$132.00
15 Foot Oxygen	CH-176	\$153.00
15 Foot Air	CH-177	\$153.00
20 Foot Oxygen	CH-178	\$174.00
20 Foot Air	CH-179	\$174.00

HUDSON RCI DISPOSABLE ADAPTERS AND CONNECTORS



No. 1667.....\$176.00
Ventilator Monitoring Adaptor Circuit

- Allows monitoring of ventilatory weaning parameters in intubated patients
- Includes oxygen nipple adaptor, aerosol tee, 6" CORR-A-FLEX tube, 15mm/22mm adaptor, m/m one-way valve and m/f one-way valve
- 50 per case



No. 1637.....\$150.00
Thermometer with Tee Adapter

- Combines No. 1635 Thermometer and No. 1077 Aerosol Tee Connector
- Packaged 50 to a case



No. 1743.....\$215.00
Valved Neb Tee

- 22mm O.D. x 18mm I.D. x 22mm I.D. connections
- Tee may be left in-line when nebulizer is removed
- Port accepts most small volume nebulizer
- 30 per case



No. 1078.....\$62.00
Ventilator Adaptor

- Right angle adapter with a capped suction port
- Vertical axis has a 15mm I.D. and horizontal axis has 22mm O.D.
- Packaged 50 to a case



No. 1641.....\$65.00
Elbow Adaptor

- Standard 22mm I.D. connector on vertical axis
- 22mm O.D. connection on horizontal axis
- Packaged 50 to a case



No. 1642.....\$84.00
Pressure Line Adaptor

- 22mm I.D./O.D. connections
- Barbed elbow accepts 5 to 7mm I.D. tubing and swivels 360 degrees
- Packaged 50 to a case



No. 1623.....\$145.00
Patient Elbow Adaptor

- 15mm O.D. x 15mm I.D./22mm O.D.
- 50 per case



No. 1648.....\$47.00
Aerosol Tee Connector

- 22mm O.D. x 15mm I.D./22mm O.D. x 22mm O.D.
- 50 per case



No. 1644.....\$107.00
One-Way Valve

- Standard 22 mm I.D./O.D. connections
- Valve offers low flow resistance and a capped monitoring port
- Packaged 50 to a case




No. 1638.....\$65.00
Aerosol Tee Connector

- Horizontal ports of the tee are 22mm O.D.
- Vertical port is 18mm I.D./22mm O.D.
- Packaged 50 to a case



No. 1077.....\$51.00
Aerosol Tee Connector

- Anti-spill design minimizes aspiration of liquids
- Horizontal axis has standard 22mm O.D. connections
- Vertical axis incorporates a standard 15mm I.D. connection and two 22mm O.D. concentric rings
- Packaged 50 to a case



No. 1639.....\$66.00
Nebulizer Tee Connector

- Standard 22 mm I.D./O.D. connections
- Nebulizer port is 18mm I.D. which accepts most small volume nebulizers
- Packaged 50 to a case



No. 1421.....\$66.00
Universal Cuff Adaptor

- Standard 22mm I.D. connections
- Packaged 50 to a case



No. 1422.....\$69.00
Multi Adaptor

- Standard 15mm I.D./22mm O.D. connections
- Packaged 50 to a case



No. 1409.....\$57.00
Large Tubing Connector

- Standard 22mm O.D. connections
- Packaged 50 to a case



No. 1640.....\$71.00
Tubing Adaptor

- Standard 22mm I.D./O.D.
- Packaged 50 to a case



No. 1751.....\$174.00
MDI Adaptor

- 15mm O.D x 15mm I.D. connections for in-line placement in 15mm tubing
- 50 per case



No. 1659.....\$167.00
MDI Adaptor

- 22mm O.D x 22mm I.D. connections for in-line placement in 22mm tubing
- 50 per case



No. 1423.....\$107.00
Tubing Adaptor

- 22mm I.D. x 5 - 7mm tapered nipple
- Packaged 50 to a case



No. 1669.....\$70.00
Monitoring Adaptor

- 15mm O.D. x 15mm I.D./22mm O.D.
- Standard luer fitting
- 50 per case

HUDSON CANNULAS & OXYGEN MASKS



- 7' over the ear
- cases of 50

CANNULAS

Part #	Description	Price
H1101	Pediatric, nonflared nasal tips (below the ears)	124.00
H1103	Adult, nonflared nasal tips	44.00
H1104	Adult, flared nasal tips	44.00



SECTION OF
STAR LUMEN TUBING

kink and crush resistant

NOSECLIP

case of 100



#64019
\$66 .00

AEROSOL MASKS



- without tubing
- cases of 50

Part #	Description	Price
H1083	Adult long	50.00
H1084	Adult short	86.00
H1085	Pedi long	70.00
H1080	Pedi short	72.00
H1095	Adult Face Tent	82.00
H1075	Adult trach	103.00
H1076	Pedi trach	125.00

NONREBREATHING



- with 7' tubing
- cases of 50

Part #	Description	Price
H1059	Adult 1 flap	121.00
H1060	Adult 2 flaps	154.00
H1058	Pedi 1 flap	166.00

MEDIUM CONCENTRATION



Part #	Description	Price
H1040	Adult short	122.00
H1041	Adult long	83.00
H1035	Pedi short	104.00
H1042	Pedi long	103.00

- with 7' tubing
- cases of 50

MULTI-VENT



- with 7' tubing
- cases of 50

- Safe, simple delivery for variable concentrations of oxygen

Part #	Description	Price
H1088	Adult	152.00
H1089	Pedi	178.00

HIGH CONCENTRATION



Part #	Description	Price
H1009	Adult short	127.00
H1007	Adult long	95.00
H1000	Pedi short	145.00
H1011	Pedi long	117.00

- with 7' tubing
- cases of 50

UNIVERSAL IPPB CIRCUIT



Universal IPPB Circuit has sensitive exhalation valve with clear cap, expiratory flow diverter, high-output 20 ml. nebulizer, mouthpiece, 6" and 39" Corr-A-Flex® tubing, multiple adaptor kit and exhalation valve drivelines

Part # 1560 \$296.00 50/cs

MDI ACCESSORY by Hudson



The POCKET CHAMBER is equipped with a highly sensitive silicone valve that opens easily without restricting the air flow. This one-way valve helps keep patients from accidentally blowing medication out of the chamber, and may aid patients experiencing difficulty timing their inhalation.

Product Description	Part #	Unit/Case	Price
Pocket Chamber Valved Aerosol Holding Chamber	1001-10	10	\$149.00

TEST LUNG (Latex Free)



Part #	Capacity	Price
LF-113	500cc	\$44.00
78175AP	1 Liter	\$44.00

AEROSOL DRAINAGE BAGS



Disposable Aerosol Drainage Bags with tee adaptor fits 22mm ID corrugated tubing. Bag has plastic chain for securing and convenient drainage port. 750 cc capacity
With wye adaptor
#1742..... \$137.00 Case of 50

DISPOSABLES



#2200 Flexible Connector (Silicone Rubber)

22mm I.D. x 22mm I.D.

This adapter was designed almost identical to the ones listed below:

INTEC 2200 MARQUEST 172200
BENNETT 3443 AM SYSTEMS 1139
Instrumentation Ind. 101-00-3443

Part Number 2200
\$247.00 per case of 50
\$7.00 each



Main Flow Bacterial / Viral Filter

- 22mm ID x 15mm ID/22mm OD
- Bacterial filtration efficiency 99.999+%
- Viral filtration efficiency 99.99+%
 - Dead space: 42mL
- Flow resistance: approximately 1.5 cmH₂O at 60 LPM
- Individually packaged 50 to a case

Part Number 1605
\$193.00 per case of 50



E-Z Wrap

Soft closed-cell foam tubes
designed to prevent ear chafing.
Light gray color
Packed 50 pair per case.

Part Number 1016
\$63.00 per case of 50



Supply Line Adapter

For extending smooth bore oxygen tubing
Tapered for excellent gripping

Part Number 65100
\$122.00 per case of 200



Corrugated Tubing

- Continuous 100-foot roll of tubing with cuttable 6-inch sections
- Packaged in protective, convenient, plastic-lined dispenser box

Part Number 1680
\$32.00 per roll



Ventilator Circuits

Adult Standard Ventilator Circuit, 72" long.
Use on Puritan-Bennett 7200A; Servo 300, 900A, 900B, and 900C; Hamilton Veolar; Emerson; Engstrom Erica; PPG Irisa; Infrasonics Adult Star and Bird 8400ST

Part Number 1613
\$94.00 per case of 10



Voldyne Incentive Spirometer

Voldyne volumetric exerciser features Good - Better - Best flow window in a compact, ergonomic design with advanced, low work-of-breathing filter design and built in handle.
5000 ml

Part Number 8884719009
\$146.00 per case of 12

NPT PIPE FITTINGS

<p>Female Tee</p>  <p>Part No. $\frac{1}{8}$" NPT..... 7034 $\frac{1}{4}$" NPT..... 7035 \$7.00</p>	<p>Street Tee</p>  <p>Part No. $\frac{1}{8}$" NPT..... 7032 $\frac{1}{4}$" NPT..... 7033 \$7.00</p>	<p>Cross</p>  <p>Part No. $\frac{1}{8}$" NPT..... 7036 $\frac{1}{4}$" NPT..... 7037 \$11.00</p>	<p>45° Street Elbow</p>  <p>Part No. $\frac{1}{8}$" NPT..... 7056 $\frac{1}{4}$" NPT..... 7057 \$7.00</p>
<p>Connectors</p>  <p>Part No. $\frac{1}{8}$" Female x $\frac{1}{8}$" Male 1017 $\frac{1}{8}$" Female x $\frac{1}{4}$" Male 7024 $\frac{1}{4}$" Female x $\frac{1}{8}$" Male..... 7025 $\frac{1}{4}$" Female x $\frac{1}{4}$" Male 7026 \$7.00</p>	<p>Hex Nipple</p>  <p>Part No. $\frac{1}{8}$" Male x $\frac{1}{8}$" Male..... 7020 $\frac{1}{4}$" Male x $\frac{1}{4}$" Male..... 7021 $\frac{1}{8}$" Male x $\frac{1}{4}$" Male..... 7022 \$7.00</p>	<p>Hose Barb Connector</p>  <p>Part No. $\frac{1}{8}$" Male x $\frac{1}{8}$" Dia. 7038 $\frac{1}{8}$" Male x $\frac{3}{16}$" Dia..... 7039 $\frac{1}{8}$" Male x $\frac{1}{4}$" Dia. 7040 $\frac{1}{4}$" Male x $\frac{1}{4}$" Dia. 7043 \$7.00</p>	<p>Coupler</p>  <p>Part No. $\frac{1}{8}$" Female x $\frac{1}{8}$" Female..... 7027 $\frac{1}{4}$" Female x $\frac{1}{4}$" Female..... 7028 $\frac{1}{4}$" Female x $\frac{1}{8}$" Female..... 7029 \$7.00</p>
<p>90° Elbow</p>  <p>Part No. $\frac{1}{8}$" Female x $\frac{1}{8}$" Female..... 7058 $\frac{1}{4}$" Female x $\frac{1}{4}$" Female..... 7059 \$7.00</p>	<p>Plug</p>  <p>Part No. $\frac{1}{4}$" NPT..... 7017 $\frac{1}{8}$" NPT..... 7019 \$7.00</p>	<p>Hose Barb Connector</p>  <p>Part No. $\frac{1}{8}$" Female x $\frac{1}{4}$" Dia. 7046 $\frac{1}{4}$" Female x $\frac{1}{4}$" Dia. 7049 \$7.00</p>	<p>Cap</p>  <p>Part No. $\frac{1}{4}$" NPT..... 7016 $\frac{1}{8}$" NPT..... 7018 \$7.00</p>
<p>Pipe Tee</p>  <p>Part No. $\frac{1}{8}$" NPT..... 7030 $\frac{1}{4}$" NPT..... 7031 \$7.00</p>	<p>Drop Ear</p>  <p>Part No. Tee $\frac{1}{8}$" F x $\frac{1}{8}$" F x $\frac{1}{8}$" F 7081 Elbow $\frac{1}{8}$" F x $\frac{1}{8}$" F..... 7082 \$18.00</p>	<p>Pipe Nipple</p>  <p>Part No. $\frac{1}{8}$" Male x $1\frac{1}{2}$" Long 7050 $\frac{1}{8}$" Male x 2" Long 7051 $\frac{1}{8}$" Male x $2\frac{1}{2}$" Long 7052 $\frac{1}{4}$" Male x $1\frac{1}{2}$" Long 7053 $\frac{1}{4}$" Male x 2" Long 7054 $\frac{1}{4}$" Male x $2\frac{1}{2}$" Long 7055 $\frac{1}{8}$" Male x $3\frac{1}{2}$" Long 7062 \$11.00</p>	<p>90° Elbow</p>  <p>Part No. $\frac{1}{8}$" F x $\frac{1}{8}$" M..... 7063 $\frac{1}{4}$" F x $\frac{1}{4}$" M..... 7065 \$7.00</p>

"Y" BLOCK ASSEMBLIES

Ohmeda Style

(Ohio)



2152	O ₂	183.00
2252	Air	183.00
2352	Vac	183.00

Chemetron Style

(NCG)



3152	O ₂	184.00
3252	Air	184.00
3352	Vac	184.00

Puritan-Bennett Style



4152	O ₂	171.00
4252	Air	171.00
4352	Vac	171.00

Ohmeda Style "Y" Assembly--"Stubby"



Compact Version

2162	O ₂	189.00
2262	Air	189.00
2362	Vac	189.00
2462	N ₂ O	189.00

T BLOCKS

OHMEDA Style



Gas	Part #	Price
Oxy	2152-T	\$197.00
Air	2252-T	\$197.00
Vac	2352-T	\$197.00



Gas	Part #	Price
Oxy	2152-TL	\$197.00
Air	2252-TL	\$197.00
Vac	2352-TL	\$197.00



Gas	Part #	Price
Oxy	2152-TR	\$197.00
Air	2252-TR	\$197.00
Vac	2352-TR	\$197.00

CHEMETRON Style



Gas	Part #	Price
Oxy	3152-T	\$197.00
Air	3252-T	\$197.00
Vac	3352-T	\$197.00



Gas	Part #	Price
Oxy	3152-TL	\$197.00
Air	3252-TL	\$197.00
Vac	3352-TL	\$197.00



Gas	Part #	Price
Oxy	3152-TR	\$197.00
Air	3252-TR	\$197.00
Vac	3352-TR	\$197.00

TRI-BLOCK ASSEMBLIES



Ohmeda Style

OR301-77-3	O ₂	257.00
OR302-77-3	Air	257.00
OR304-77-3	Vac	257.00



Chemetron Style

OR301-76-3	O ₂	239.00
OR302-76-3	Air	239.00
OR304-76-3	Vac	239.00

DISS "Y" ASSEMBLY WITH CHECK VALVES



(M24-77)



(M24-76)



(7211)

7211	Diss Female	O ₂	86.00
M24-76	Chemetron	O ₂	94.00
M24-77	Ohmeda (Ohio)	O ₂	95.00
M24-79	Puritan	O ₂	79.00
7212	Diss Female	Air	86.00
M16-76	Chemetron	Air	101.00
M16-77	Ohmeda (Ohio)	Air	101.00
M16-79	Puritan	Air	89.00
7213	Diss Female	Vac	86.00
M22-76	Chemetron	Vac	101.00
M22-77	Ohmeda (Ohio)	Vac	101.00
M22-79	Puritan	Vac	89.00

FOR MORE "Y" BLOCK ASSEMBLIES

PLEASE VISIT OUR WEBSITE

WWW.RCMEDICAL.COM

QUICK CONNECTS (MALE)

Ohmeda

1/8" NPT Male



Gas	Part #	Price
Oxygen.....	2101.....	\$39.00
Air.....	2201.....	39.00
Vacuum.....	2301.....	39.00

1/8" NPT Female



Oxygen.....	2102.....	\$37.00
Air.....	2202.....	37.00
Vacuum.....	2302.....	37.00

1/4" NPT Male



Oxygen.....	2103.....	\$38.00
Air.....	2203.....	38.00
Vacuum.....	2303.....	38.00

DISS Male



Oxygen.....	2105.....	\$45.00
Air.....	2205.....	45.00
Vacuum.....	2305.....	45.00

DISS Hex Nut



Oxygen.....	2106.....	\$45.00
Air.....	2206.....	45.00
Vacuum.....	2306.....	45.00

DISS Hand Tight



Oxygen.....	2108.....	\$46.00
Air.....	2208.....	46.00
Vacuum.....	2308.....	46.00

1/4" Hose Barb



Oxygen.....	2111.....	\$44.00
Air.....	2211.....	44.00
Vacuum.....	2311.....	44.00

DISS Male with Check Valve



Oxygen.....	2114.....	\$48.00
Air.....	2214.....	48.00
Vacuum.....	2314.....	48.00

Chemetron

1/8" NPT Male



Gas	Part #	Price
Oxygen.....	3101.....	\$35.00
Air.....	3201.....	35.00
Vacuum.....	3301.....	35.00

1/8" NPT Female



Oxygen.....	3102.....	\$35.00
Air.....	3202.....	35.00
Vacuum.....	3302.....	35.00

1/4" NPT Male



Oxygen.....	3103.....	\$37.00
Air.....	3203.....	37.00
Vacuum.....	3303.....	37.00

DISS Male



Oxygen.....	3105.....	\$42.00
Air.....	3205.....	42.00
Vacuum.....	3305.....	42.00

DISS Hex Nut



Oxygen.....	3106.....	\$42.00
Air.....	3206.....	42.00
Vacuum.....	3306.....	42.00

DISS Hand Tight



Oxygen.....	3108.....	\$44.00
Air.....	3208.....	44.00
Vacuum.....	3308.....	44.00

1/4" Hose Barb



Oxygen.....	3111.....	\$39.00
Air.....	3211.....	39.00
Vacuum.....	3311.....	39.00

DISS Male with Check Valve



Oxygen.....	3114.....	\$45.00
Air.....	3214.....	45.00
Vacuum.....	3314.....	45.00

Puritan Bennett

1/8" NPT Male



Gas	Part #	Price
Oxygen.....	4101.....	\$33.00
Air.....	4201.....	33.00
Vacuum.....	4301.....	33.00

1/8" NPT Female



Oxygen.....	4102.....	\$33.00
Air.....	4202.....	33.00
Vacuum.....	4302.....	33.00

1/4" NPT Male



Oxygen.....	4103.....	\$35.00
Air.....	4203.....	35.00
Vacuum.....	4303.....	35.00

DISS Male with Knob



Oxygen.....	4105.....	\$40.00
Air.....	4205.....	40.00
Vacuum.....	4305.....	40.00

DISS Hex Nut



Gas	Part #	Price
Oxygen.....	4106.....	\$37.00
Air.....	4206.....	37.00
Vacuum.....	4306.....	37.00

DISS Hand Tight



Oxygen.....	4108.....	\$39.00
Air.....	4208.....	39.00
Vacuum.....	4308.....	39.00

1/4" Hose Barb w/ Knob



Oxygen.....	4111.....	\$39.00
Air.....	4211.....	39.00
Vacuum.....	4311.....	39.00

DISS Male w/ Knob & Check Valve



Oxygen.....	4114.....	\$44.00
Air.....	4214.....	44.00
Vacuum.....	4314.....	44.00

QUICK CONNECTS (FEMALE)

CHEMETRON STYLE CHECK UNITS (Female Couplers)

1/8" NPT Female



Gas	Part No.
O ₂	OR 101
Air.....	OR 102
Vac*.....	OR 104

DISS Male



Gas	Part No.
O ₂	OR 101-2
Air.....	OR 102-2
Vac*.....	OR 104-2

1/4" NPT Male



Gas	Part No.
O ₂	OR 101-3
Air.....	OR 102-3
Vac*.....	OR 104-3

DISS Female



Gas	Part No.
Hand Tight O ₂	OR 101-5
Hand Tight Air.....	OR 102-5
Hand Tight Vac*.....	OR 104-5

\$72.00

OHMEDA STYLE CHECK UNITS (Female Couplers)

1/8" NPT Female



Gas	Part No.
O ₂	OR 201
Air.....	OR 202
Vac*.....	OR 204

DISS Male



Gas	Part No.
O ₂	OR 201-2
Air.....	OR 202-2
Vac*.....	OR 204-2

1/4" NPT Male



Gas	Part No.
O ₂	OR 201-3
Air.....	OR 202-3
Vac*.....	OR 204-3

DISS Female



Gas	Part No.
Hand Tight O ₂	OR 201-5
Hand Tight Air.....	OR 202-5
Hand Tight Vac*.....	OR 204-5

\$78.00

PURITAN-BENNETT STYLE CHECK UNITS (Female Couplers)

1/8" NPT Female



Gas	Part No.
O ₂	4150
Air.....	4250
Vac*.....	4350

DISS Male



Gas	Part No.
O ₂	4155
Air.....	4255
Vac*.....	4355

1/4" NPT Male



Gas	Part No.
O ₂	4153
Air.....	4253
Vac*.....	4353

DISS Female



Gas	Part No.
Hand Tight O ₂	4158
Hand Tight Air.....	4258
Hand Tight Vac*.....	4358

\$69.00

*Vacuum primary check only

SCHRADER STYLE CHECK UNITS (Female Couplers)

1/8" NPT Female



Gas	Part No.
O ₂	6150
Air.....	6250
Vac*.....	6350

DISS Male



Gas	Part No.
O ₂	6155
Air.....	6255
Vac*.....	6355

1/4" NPT Male



Gas	Part No.
O ₂	6153
Air.....	6253
Vac*.....	6353

DISS Female



Gas	Part No.
Hand Tight O ₂	6158
Hand Tight Air.....	6258
Hand Tight Vac*.....	6358

\$94.00

DISS FITTINGS

DISS 90° Elbow Male x Female



Oxygen	Air	Vac
#0705	#0705-A	#0705-V
\$16.00	\$42.00	\$42.00

DISS Male x NPT Male



Oxygen	Air	Vac
1/8" #0701...\$6.00	#0731...\$15.00	#0751...\$15.00
1/4" #M24-5...\$8.00	#0733...\$15.00	#0753...\$15.00

DISS Male x DISS Male



Oxygen	Air	Vac
#0704	#0704-A	#0704-V
\$8.00	\$42.00	\$42.00

DISS Male x 1/8" NPT Male with Check Valve



Oxygen	Air	Vac
1/8" #0712...\$13.00	#0742...\$14.00	#0762...\$22.00
1/4" #0713...\$13.00	#0743...\$14.00	#0763...\$22.00

DISS Female x DISS Female



Oxygen	Air	Vac
#0104	#0104-A	#0104-V
\$17.00	\$29.00	\$29.00

DISS Male x DISS Female



Oxygen	Air	Vac
#1005-SP	#1005-SPA	#1005-SPV
\$35.00	\$35.00	\$35.00

DISS Female x 1/8" NPT Male



Oxygen	Air	Vac
#0110	#0209	#0309
\$10.00	\$13.00	\$15.00

DISS Male x 1/8" NPT Female



Oxygen	Air	Vac
#0702	#0702A	#0702V
\$6.00	\$16.00	\$16.00





50 PSI
Auxiliary
Valve

OPTO - \$16.00
(OXYGEN)

APTO - \$19.00
(AIR)


Thread Comparisons Guide:



(Measurements given are approximate.)

For Oxequip, Medstar and Schrader
quick connects

PLEASE VISIT OUR WEBSITE AT:
WWW.RCMEDICAL.COM



NEB DRIVE RESCUE FLOW

An Economical adapter that delivers a preset flow from standard wall outlets.
It's a convenient, dependable way to administer aerosols.

DISS FEMALE



\$30.00

O.	6 to 7 LPM	8116
O.	7 to 8 LPM	8117
O.	8 to 9 LPM	8118
O.	10 to 11 LPM	8119
Air	6 to 7 LPM	8216
Air	7 to 8 LPM	8217
Air	8 to 9 LPM	8218
Air	10 to 11 LPM	8219
Rescue Flow	14 to 16 LPM	8115
O.	24 to 26 LPM	8114

OHMEDA



\$49.00

O.	6 to 7 LPM	8126
O.	7 to 8 LPM	8127
O.	8 to 9 LPM	8128
O.	10 to 11 LPM	8129
Air	6 to 7 LPM	8226
Air	7 to 8 LPM	8227
Air	8 to 9 LPM	8228
Air	10 to 11 LPM	8229
Rescue Flow	14 to 16 LPM	8125
O.	24 to 26 LPM	8124

CHEMETRON



\$49.00

O.	6 to 7 LPM	8136
O.	7 to 8 LPM	8137
O.	8 to 9 LPM	8138
O.	10 to 11 LPM	8139
Air	6 to 7 LPM	8236
Air	7 to 8 LPM	8237
Air	8 to 9 LPM	8238
Air	10 to 11 LPM	8239
Rescue Flow	14 to 16 LPM	8135
O.	24 to 26 LPM	8134


PURITAN-BENNETT



\$35.00

O.	6 to 7 LPM	8146
O.	7 to 8 LPM	8147
O.	8 to 9 LPM	8148
O.	10 to 11 LPM	8149
Air	6 to 7 LPM	8246
Air	7 to 8 LPM	8247
Air	8 to 9 LPM	8248
Air	10 to 11 LPM	8249
Rescue Flow	14 to 16 LPM	8145
O.	24 to 26 LPM	8144

To order Neb Drive or Rescue Flow with
shut off valve, please add an "S" to the end
of the part number.
The price is \$40.00 more.



Shut off valve allows you to leave neb drive
in wall outlet.

REPLACEMENT RESPIRATORY OXYGEN SENSORS



Part # SS-05

OEM Model / Sensor:
Hudson 5577 / 5500
Hudson 6477 / 5500
Hudson 5590 / 5500
Teledyne R13



Part # SS-03

OEM Model / Sensor:
Newport OM100 / T-7
Newport VM200 / T-7
Sechrist 650 / 43170
Teledyne 60T / T-7
Teledyne 191 / T-7
Teledyne 200 / T-7
Hudson / 5558

Fits ALL MiniOx!



Part # SS-11A

OEM Model / Sensor:
Bird 6400 / 66030
Bird 8400 / 66030
Lifecare PLV-102 / 37-079
All MiniOx/406931
Puritan Bennett 7200 Vent
& 7820 Monitor / 4-020933
Respironics Esprit / 1001454
Hudson / 5567



Part # SS-25A

OEM Model / Sensor:
Maxtec OM-25A / 250E
Maxtec OM-25ME / 250E



Part # SS-12A

OEM Model / Sensor:
Hamilton Aladdin / 51254
Hamilton Galileo / 51254
Hamilton Raphael / 51254
Infrasonic Adult Star / 6011
Marquet Saracap
Puritan Bennett 760 & 840
Ventilators / 4-074925-00
Siemens 300 / 6419332
Siemens 900 / 9004979
Siemens Servo 710 / 9004979
Spacelab / 90514
Hudson / 5566
Replaces
Teledyne R22
MSA 472062


ALL SENSORS \$127.00

Unconditional 12 month Warranty!
Exceptional Quality and Performance!
Canned for increased shelf life and freshness!




**FOR MORE SENSORS
PLEASE VISIT OUR WEBSITE
WWW.RCMEDICAL.COM**

MINIOX 200E OXYGEN ANALYZER



The MiniOX 200E is best used anywhere an inline analyzer is needed, when alarms are not deemed important. It is especially useful in conjunction with older ventilators or oscillators. Designed with simplicity and portability in mind, this analyzer quickly adapts for use in tubing applications and is ideal for checking oxygen percent accuracy in line or on concentrations. Its patented ergonomic design fits comfortably in your hand and features one-touch on/off and calibration. The MiniOX 200E will calibrate to room air (20.9%) or 100% Oxygen. It is ideal for applications needing an external sensor.

- 1 Large LCD - Easy-to-read characters
- 2 Calibration Reminder - Reminds once a week
- 3 One-Touch Calibration - Auto calibrate to room air or 100% Oxygen
- 4 On/Off Button - One-touch on and off
- 5 Patented Ergonomic Design - Fits comfortably in your hand
- 6 External Sensor Cable - 10 ft. extended, 3 ft. retracted
- 7 Low Battery Indicator - 5,000+ Hours of battery life
- 8 Rubber Gasket - Prevents moisture from contact board
- 9 External Sensor - MiniOX 250E, 0-100% Oxygen, 24 month warranty



MiniOX 200E Oxygen Analyzer

Part # 600210..... \$629.00

Comes complete with sensor, cable, tee and AA batteries

Replacement Sensor # 606940 \$179.00

Approximately 5,000 hours with continuous use

MISCELLANEOUS

PLASTIC BAGS (DISPOSABLE)



PATIENT SET-UP (Cotton Drawstring)

Part No.	Description	Quantity	Price
0049	8 x 10	500	\$109.00
0059	11 3/4 x 15	500	\$117.00
XL-S2	16 x 18 x 3	250	\$112.00

EQUIPMENT COVERS

Dust covers designed to protect ventilators and other equipment. Comes on a roll.

Part No.	Description	Quantity	Price
PLC-60	28 x 22 x 58 (clear)	55	\$109.00
PLC-60/B	28 x 22 x 58 (buff)	55	\$116.00
0084	28 x 22 x 56 (clear***)	50	\$105.00
7083	28 x 22 x 56 (yellow)	50	\$105.00
0126	48 x 25 x 42 (red)	50	\$124.00

***printed, "CLEAN, READY FOR USE"

POSEY CUFFLATOR™

ENDOTRACHEAL TUBE INFLATOR & MANOMETER



An accurate, convenient tracheal cuff inflator and manometer. No more syringes, electricity or complicated connections needed to check cuff pressure regularly. Air vent button and inflator bulb let you vary cuff pressure quickly and easily. Clip on back allows easy positioning. Gauge shows recommended pressure range in cm H₂O.

#8199 Cufflator with instrument case.....\$469.00

RUSCH POSITIVE - NEGATIVE PRESSURE GAUGE

60 - 0 - 60 mm. hg.

Measures: • Millimeters of Mercury
• Centimeters of Water



(395943)

33-3505



Hose Barb Fitting

395943 Gauge only\$85.00

33-3505 Tubing, T piece,\$10.00 each or
Trach fitting \$329.00/case of 48

ANESTHESIA ASSOCIATES PRESSURE MANOMETER



274

Part#	Description	Price
199	Gauge with 1/4" NPT	\$91.00
274	Gauge with 5-55 cm H ₂ O adjustable pressure relief valve and "T" for supply tubing.	152.00

Part#	Description	Price
199-MAX	Gauge with Pointer	139.00
274-MAX	274 with Pointer	179.00


Manometer has range of -40 to +80 cm H₂O in a 2 1/2" diameter chrome plated case with 1/4" NPT bottom connection. Has front zero correcting screw.



199-MAX

MAX is a resettable low friction pointer to capture maximum pressure or vacuum reading

PLEASE VISIT OUR WEBSITE
WWW.RCMEDICAL.COM




Peak Flowmeter

Part number232-6000
Case of 10\$227.00



Range: 60 - 850 LPM
Scale Intervals: 10 LPM
Reading Accuracy: 5 LPM
Measurement Accuracy: +/- 10%
Repeatability: +/- 5%



Hinged handle for patient convenience

COMPRESSORS AND NEBULIZERS



Uniquely designed for bed-side use. Special material absorbs most noise for quiet performance. Audible and visual high temperature and low pressure alarms with DC powered six-volt backup. Vibration-free gauge shows diminished pressure, impaired drying or exceeded flow capacity. Circuit breaker and thermal overload protector included. Continuous flow of 64 LPM at 50 psi.

Part Number TE20 \$4550.00



Designed for bedside use. Quiet, vibration-free performance. Provides condensation-free compressed air to respiratory therapy equipment. Vibration-free gauge shows diminished pressure, impaired drying or exceeded flow capacity. Circuit breaker and thermal overload protector included. Continuous flow of 42 LPM at 50 psi.

Part Number TE35 \$3880.00



DeVilbiss Pulmo-Aide® Medication Nebulizer

Durable AC-powered medication nebulizer with Pulmo-Neb® 60-day disposable nebulizer and tubing. Unit weighs 7 lbs. and has 9 - 16 LPM range. 5 year warranty. 115 vac

3-Prong Plug

DV4650D \$180.00

BLENDERS - PRECISION MEDICAL

Air/O₂



Its unique. It is the only unit available that allows the clinician to turn the bleed control to ON or OFF providing quieter operation. This feature is especially useful in the NICU and PICU.

Part #	Flow	Rate	Price
5200	High Flow	15 - 120 lpm	999.00
5300	Low Flow	3 - 30 lpm	999.00

Heliox



Its unique. It is the only free-standing pneumatic HeliO₂ blending device on the market today. The ability to turn the bleed control collar ON or OFF provides ± 3% F₁₀₂ accuracy as well as quieter operation and maximal Heliox tank durations.

Part #	Flow	%	Price
5570	High Flow	70 - 30	1795.00
5580	High Flow	80 - 20	1795.00
5470	Low Flow	70 - 30	1795.00
5480	Low Flow	80 - 20	1795.00

ARTERIAL BLOOD GAS KITS

Pulset cricket system with needle, 25U balance heparin. Alcohol pad, iodophor preps, 2 x 2 gauze, ice bag, bandage and patient ID.



Part # 336295
Kit 22G x 1" N
\$490.00/Case 100

Part # 336395
Kit 23G x 1" N
\$432.00/Case 100

RC MEDICAL SELLS A WIDE VARIETY OF RESPIRATORY SUPPLIES AND EQUIPMENT. PLEASE VISIT OUR WEBSITE AT: WWW.RCMEDICAL.COM

REGULATOR AND FLOWMETER ACCESSORIES



7-5A



7-5C



0159



0159C

Nut and Nipple (Christmas Tree)

Part #	Description	Price
7-4A*	One piece Green	\$27.00 cs/50
7-5A*	Two piece Green	\$34.00 cs/50
7-5C*	Two piece Yellow	\$46.00 cs/50
7-5H*	Two piece Clear	\$46.00 cs/50

Plastic - For single Patient Use (Above)*

0159	Metal Nut/Nipple	\$ 7.00 each
0159C	Metal w/ 6" chain	\$ 8.00 each

CYLINDER WRENCHES FOR SMALL AND LARGE CYLINDERS



792-A



792-B



792-D



MCW-2B



66071



MCW-2BC



WR-3



5080

Part #	Description	Price
7-92A	Plastic Small Wrench, BLACK	26.00 cs/10
7-92B	Plastic Small Wrench, GREEN	26.00 cs/10
7-92D	Plastic Small Wrench, PINK	26.00 cs/10
66071	Plastic Small Hand Wheel	65.00 cs/10
MCW-2B	Metal Small Cylinder Wrench	36.00 cs/10
MCW-2BC	Metal Small Wrench with Chain	46.00 cs/10
5080	Hudson Small Metal Wrench	8.00 each
WR-3	Metal Wrench Small & Large	12.00 each

REPLACEMENT GAUGES FOR OXYGEN REGULATORS

2" Bottom Port 1/4" NPT



2510




2530



2543

1 1/2" Back Port 1/8" NPT



2501



2521

Low Profile 3000 PSI



MJ-C1

Part #	Description	Price
2510	15 lpm - 2" Bottom Port	\$25.00
2530	Contents (0-4000 PSI) 2" Bottom	\$25.00
2543	100 PSI 2" Bottom Port	\$25.00
2501	15 lpm - 1 1/2" Back Port	\$25.00
2521	Contents (0-4000 PSI) 1 1/2" Back	\$32.00
2500	15 lpm - 1 1/2" Bottom Port	\$25.00
2520	Contents (0-4000 PSI) 1 1/2" Bottom	\$25.00
MJ-C1	Low Profile Contents (3000 psi)	\$15.00
AIR, CO2 GAUGES		
MG-C11	AIR, CO2 (0 - 100 PSI) 2" Bottom	\$39.00
MG-C9	AIR, CO2 (0-4000 PSI) 2" Bottom	\$39.00

Yoke Washer For Small Cylinder Regulators (CGA 870)



2544-2

Part #	Description	Price
2544-2	aluminum w/ rubber center	\$ 3.00 each
2544-2C	aluminum w/ rubber center	\$98.00 cs/50

LITERMETER

Every therapist should have one for checking the accuracy of regulators and flowmeters.



Part #	LPM	Price
1271	0-2.5	\$37.00
MT200	0-8	37.00
MT300	2-15	37.00

OXYGEN TIMER



TIMOX features:

- Automatic timing of oxygen flow as low as .30 lpm.
- Unit easily reset to zero with tip of pen or pencil.
- 1999.9 recorded hours of oxygen flow.
- Two year warranty on unit excluding battery.

Timer Part #8050 \$133.00
Replacement Battery Part #8066 \$18.00

LO₂ PRESSURE GAUGE (Low Pressure Alarm)

**Audible and Visible Alarm
When Tank Pressure Falls Below 500 PSI**



Part #	Description	Price
S2000-1/4	1/4" Thread	\$276.00
S2000-1/8	1/8" Thread	\$276.00

RC Medical, Inc.
284 Merrow Road
P.O. Box 833
Tolland, CT 06084

Please deliver this catalog to the Purchasing
Department if addressee is not available

ALL IN STOCK ORDERS PLACED BY 2 PM (EST) SHIPPED THE SAME DAY.



Ohmeda 50 PSI Oxygen
Regulator with
Chemetron Female
Coupler
LQ-1..... \$189.00



Ohmeda 50 PSI Oxygen
Regulator with
Ohmeda Female Coupler
LQ-2..... \$189.00



Ohmeda 50 PSI Oxygen
Regulator with Puritan
Bennett Female Coupler
LQ-3..... \$189.00



Ohmeda 50 PSI Oxygen
Regulator with
DISS Male Outlet
LQ-4..... \$129.00

PLEASE VISIT OUR WEBSITE AT: WWW.RCMEDICAL.COM



Ohmeda Medical 15
LPM Oxygen Regulator
for Small Cylinder
(CGA 870)
LQ-5..... \$149.00




Victor 50 PSI H Tank
Regulator
LQ-6..... \$59.00




Allied Healthcare 15 LPM
Oxygen Regulator for
Large Cylinder (CGA 540)
LQ-7..... \$99.00



Chemetron 0 - 15 LPM
Oxygen Regulator, H
Cylinder
LQ-8..... \$79.00



FREE GAUGE SAVERS
ON ABOVE REGULATORS
MADE IN THE USA
LIMITED QUANTITIES



PHONE: (888) 431-2727

FAX: (888) 337-0443