RESPIRATORY CATALOG

EVERYDAY AND SPECIALTY PRODUCTS







Order online at: www.rcmedical.com





RC MEDICAL, INC.

Phone: 888-431-2727 • Fax: 888-337-0443

2023

P.O. Box 833

284 Merrow Road

Tolland, CT 06084

HEALTHCARE FACILITIES HAVE AN OPEN 30 DAY ACCOUNT

NO MINIMUM ORDER - PURCHASE ORDERS ACCEPTED

F.O.B. Origin. • Prices subject to change without notice.

25% restocking charge prior to 30 days. • Please call for a RGA Number.

Products may NOT be returned after 30 days.

PRICING ON WEBSITE SUPERSEDES CATALOG PRICES

INDEX

ADAPTERS	
	HOSES17, 18
DISS27	MISCELLANEOUS29, 30
Female Couplers25	OXYGEN ANALYZERS AND SENSORS28
Male Quick Connects24	RAIL SYSTEMS16
NPT Fittings22	REGULATORS
Y Blocks and T Blocks23	Air15
BAGS, EQUIPMENT COVERS29	CO215
COMPRESSORS30	Heliox15
CYLINDER CARTS AND STANDS9 - 11	Oxygen12, 13
FITTINGS (under ADAPTERS)	Oxygen - Specialty14
FLOWMETERS5 - 8	Suction2 -4
FLOWMETER ACCESSORIES31	RESPIROMETERS21

SUCTION REGULATOR SPECIALS LIMITED QUANTITY



	Ohmeda Wall Adapter	Chemetron Wall Adapter	Puritan Wall Adapter	DISS Wall Adapter
		==		;
	Part #	Part #	Part #	Part #
Ohio Continuous	LQ-49	LQ-51	LQ-52	LQ-50
Allied Intermittent	LQ-80	LQ-81	LQ-82	LQ-83
Amvex Intermittent	LQ-21	LQ-22	LQ-23	LQ-24

All suction regulators above come standard with tubing nipple outlet.



For DISS Male outlet, add a "D" to end of part number.

Toll Free Phone 888-431-2727 1-860-871-4156

Toll Free Fax 888-337-0443 1-860-870-0443







SUCTION REGULATORS – by Precision Medical CONTINUOUS AND CONTINUOUS/INTERMITTENT

W LOWER PRICING

MADE IN USA 10 YEAR WARRANTY **USER FRIENDLY LARGE KNOB** EASY TO READ FACE

2 Mode Continuous \$279





Wall Adapter Inlet **DISS Female** Ohmeda Chemetron

Puritan

Tubing Nipple Outlet PM3003 PM3004 PM3005 Tubing PM3008 Nipple

3 Mode Continuous/Full line \$279





Wall Adapter Inlet **DISS Female** Ohmeda Chemetron **Puritan**

Ohmeda

Puritan

Chemetron







3 Mode Continuous/Intermittent \$398 Wall Adapter **Tubing Nipple** Inlet Outlet **DISS Female** PM3303



PM3304 PM3305 PM3308 Nipple

DISS Male Outlet PM3314 PM3315 PM3316 DISS Male PM3319

WALL ADAPTERS (INLET)



DISS Female w/hand nut



Ohio (Ohmeda



Chemetron (NCG)



Puritan

VACUTRON SUCTION REGULATORS by Allied

The new, sleekly contoured Vacutron[©] is a precision engineered instrument that combines the most desirable features for suction regulators. Smaller in size, but big on benefits, it offers advantages that include:





Easy access ON/OFF Soft-Touch knob gives full-line capability on the surgical model.

Maintenance is easy. Clean by flushing with cold sterilant.

Weight-bearing brass post prevents case from cracking. Various inlet and outlet fittings insure compatibility with existing fixtures in any

for durability.

0-200 mmHg is achieved easily with minimal turning of the knob. Vacuum-limiting valve on the continuous model limits vacuum at 300 mmHg to prevent potential tissue damage.

New larger Soft-Touch knob for quick, fumble-free action. Exposed bottom allows thumb-wheel adjustment. Oversized knob for easy operation.

Price

\$499.00

499.00

499.00

499.00



Part #

22-13-1106

22-13-1107

22-13-1108

22-13-1105

22-13-1206

22-13-1207

22-13-1208

22-13-1205

medical facility.

Diss Female w/hand nut



Chemetron (NCG)



Ohio (Ohmeda)



Puritan



Tubing Nipple



Continuous/Intermittent Suction Regulator			
22-15-1106	Ohmeda	DISS Male	\$899.00
22-15-1107	DISS Female	DISS Male	899.00
22-15-1108	Chemetron (NCG)	DISS Male	899.00
22-15-1105	Puritan-Bennett	DISS Male	899.00

Continuous Suction Regulator

Outlet

DISS Male

DISS Male

DISS Male

DISS Male

Tubing Nipple 499.00

Tubing Nipple 499.00 Tubing Nipple 499.00

Tubing Nipple 499.00

Inlet

Ohmeda

Ohmeda

DISS Female

Chemetron (NCG)

Puritan-Bennett


DISS Female

Chemetron (NCG)

Puritan-Bennett

22-15-1206	Ohmeda	Tubing Nipple	899.00
22-15-1207	DISS Female	Tubing Nipple	899.00
22-15-1208	Chemetron (NCG)	Tubing Nipple	899.00
22-15-1205	Puritan-Bennett	Tubing Nipple	899.00

Collection Canisters



Disposable Collection Canisters

Size	Outlet	Units/Case	Part Number	Price
2400ML	DISS	36/CASE	20-08-0001	\$294.00
2400ML	STEM	36/CASE	20-08-0002	287.00
1500ML	STEM	48/CASE	20-08-0003	306.00
1500ML	DISS	48/CASE	20-08-0004	313.00
2400ML	DISS	36/CASE w/tubing	20-08-0005	368.00
2400ML	STEM	36/CASE w/tubing	20-08-0006	387.00
1500ML	STEM	42/CASE w/tubing	20-08-0007	422.00
1500ML	DISS	42/CASE w/tubing	20-08-0008	409.00

Suction Canister Holder and Wall Bracket





White Plastic Holder 20-02-0150 \$34.00

Wall Bracket 22-11-0006 \$72.00

For suction accessories please visit our website at:

www.rcmedical.com

Tubing Fittings Vac Trap

Stands Bacteria Filters Wye Blocks

PRECISION MEDICAL CHROME BODY FLOWMETER

5 Year Warranty



4MFA1001



1MFA9001



3MFA1001



0 - 1 LPM PEDIATRIC OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
4MFA1001	None - 1/8" NPT Female	\$239.00
4MFA1003	DISS FEMALE w/Hand Tight	252.00
4MFA1006	CHEMETRON (NCG)	257.00
4MFA1005	OHIO/OHMEDA	261.00
4MFA1008	PURITAN BENNETT	252.00

0 - 3 LPM PEDIATRIC OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA3001	None - 1/8" NPT Female	\$103.00
1MFA3003	DISS FEMALE w/Hand Tight	110.00
1MFA3006	CHEMETRON (NCG)	134.00
1MFA3005	OHIO/OHMEDA	134.00
1MFA3008	PURITAN BENNETT	112.00

0 - 15 LPM OXYGEN FLOWMETER WITH ADAPTER

PART NU	MBER	ADAPTER	PRICE
1MFA1	001	None - 1/8" NPT Female	\$70.00
1MFA1	003	DISS FEMALE w/Hand Tight	79.00
1MFA1	006	CHEMETRON (NCG)	99.00
1MFA1	005	OHIO/OHMEDA	99.00
1MFA1	008	PURITAN BENNETT	90.00

0 - 15 LPM OXYGEN FLOWMETER WITH 50 PSI AUXILIARY VALVE AND ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA1003PTO	DISS FEMALE w/Hand Tight	\$106.00
1MFA1006PTO	CHEMETRON (NCG)	131.00
1MFA1005PTO	OHIO/OHMEDA	131.00
1MFA1008PTO	PURITAN BENNETT	114.00
1MFA100MPTO	MEDSTAR	128.00

0 - 70 LPM HI FLOW OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA8001	None - 1/8" NPT Female	\$113.00
1MFA8003	DISS FEMALE w/Hand Tight	119.00
1MFA8006	CHEMETRON (NCG)	131.00
1MFA8005	OHIO/OHMEDA	132.00
1MFA8008	PURITAN BENNETT	124.00

0 - 70 LPM HI FLOW AIR FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA9001	None - 1/8" NPT Female	\$113.00
1MFA9003	DISS FEMALE w/Hand Tight	126.00
1MFA9006	CHEMETRON (NCG)	130.00
1MFA9005	OHIO/OHMEDA	133.00
1MFA9008	PURITAN BENNETT	127.00

SELECT FLOWMETER

15 LPM OXYGEN FLOWMETER WITH ADAPTER WITH FLOW SELECTOR BUILT IN

PART NUMBER	ADAPTER	PRICE
3MFA1001	None - 1/8" NPT Female	\$87.00
3MFA1003	DISS FEMALE w/Hand Tight	95.00
3MFA1006	CHEMETRON (NCG)	105.00
3MFA1005	OHIO/OHMEDA	107.00
3MFA1008	PURITAN BENNETT	99.00

FOR MORE FLOWMETERS, PLEASE VISIT OUR WEBSITE AT: WWW.RCMEDICAL.COM


FLOW SELECTOR



Flow Selector Features

- Give nebulizer treatments without disconnecting the humidifier
- Reduces the changeover time to seconds
- No more loose tubing laying around
- Eliminates tubing nipple adaptors (It has 2 of its own.)
- Fits all flowmeters that have DISS male outlet fittings
- Lightweight

Part # Price 1000 \$52.00



O-200cc, 0-1LPM, 0-3 LPM, 0-70 LPM, 0-15 LPM

5

PRECISION MEDICAL COMPACT 15 LPM FLOWMETERS

8MFA1005



8MFA1006PTO



Y8MFA1006



Y8MFA1005PT2

OXYGEN

15 LPM OXYGEN FLOWMETER WITH ADAPTER

Part #	Adapter	Price
8MFA1001	None - 1/8" NPT Female	\$45.00
8MFA1002	Diss Female	51.00
8MFA1003	Diss Female w/hand nut	52.00
8MFA1006	Chemetron (NCG)	76.00
8MFA1005	Ohio (Ohmeda)	78.00
8MFA1008	Puritan	64.00
8MFA100M	Medstar	103.00
8MFA1007	Oxequip	61.00

15 LPM OXYGEN FLOWMETER WITH 50 PSI AUXILIARY VALVE AND ADAPTER

Part #	Adapter	Price
8MFA1001 PT0	None - 1/8" NPT Female	\$60.00
8MFA1002 PT0	Diss Female	62.00
8MFA1003 PT0	Diss Female w/hand nut	79.00
8MFA1006 PT0	Chemetron (NCG)	85.00
8MFA1005 PT0	Ohio (Ohmeda)	85.00
8MFA1008 PT0	Puritan	76.00
8MFA100M PT0	Medstar	112.00
8MFA1007 PT0	Oxequip	73.00

15 LPM DUAL OXYGEN FLOWMETER WITH ADAPTER

Part #	Adapter	Price
Y8MFA1001	None - 1/8" NPT Female	\$97.00
Y8MFA1002	Diss Female	102.00
Y8MFA1003	Diss Female w/hand nut	104.00
Y8MFA1006	Chemetron (NCG)	127.00
Y8MFA1005	Ohio (Ohmeda)	126.00
Y8MFA1008	Puritan	111.00
Y8MFA100M	Medstar	154.00
Y8MFA1007	Oxequip	112.00

15 LPM DUAL OXYGEN FLOWMETER WITH (2) 50 PSI AUXILIARY VALVES AND ADAPTER

Part #	Adapter	Price
Y8MFA1001 PT2	None - 1/8" NPT Female	\$ 109.00
Y8MFA1002 PT2	Diss Female	114.00
Y8MFA1003 PT2	Diss Female w/hand nut	115.00
Y8MFA1006 PT2	Chemetron (NCG)	146.00
Y8MFA1005 PT2	Ohio (Ohmeda)	144.00
Y8MFA1008 PT2	Puritan	122.00
Y8MFA100M PT2	Medstar	165.00
Y8MFA1007 PT2	Oxequip	123.00

CHAIN-RING ASSEMBLY



The Chain-Ring Assembly is used to permanently attach Flowmeters to the wall outlet.

Part #	Description	Price
8004	6" Chain w/Closed Ring	\$7.00
8005	6" Chain w/Open Ring	\$7.00



Diss Female



Diss Female with hand nut



Chemetron (NCG)



Ohio (Ohmeda)



Puritan



Medstar



Oxequip



50 PSI Auxiliary Valve

PRECISION MEDICAL COMPACT 15 LPM FLOWMETERS

AIR

Part #

8MFA2001 PT0

8MFA2002 PT0

8MFA2003 PT0

8MFA2006 PT0

8MFA2005 PT0

8MFA2008 PT0

8MFA200M PT0

8MFA2007 PT0



8MFA2005

15 LPM AIR FLOWMETER WITH ADAPTER		
Part #	Adapter	Price
8MFA2001	None - 1/8" NPT Female	\$45.00
8MFA2002	Diss Female	54.00
8MFA2003	Diss Female w/hand nut	59.00
8MFA2006	Chemetron (NCG)	79.00
8MFA2005	Ohio (Ohmeda)	89.00
8MFA2008	Puritan	68.00
8MFA200M	Medstar	107.00
8MFA2007	Oxequip	63.00

15 LPM AIR FLOWMETER

WITH 50 PSI AUXILIARY VALVE AND ADAPTER

None - 1/8" NPT Female

Diss Female w/hand nut

Chemetron (NCG)

Ohio (Ohmeda)

Adapter

Puritan

Medstar

Oxequip

Diss Female





Diss Female with hand nut

Price

\$61.00

67.00

74.00

89.00

99.00

84.00

122.00

79.00



Chemetron (NCG)



Ohio (Ohmeda)



Puritan



Medstar



Oxequip




50 PSI Auxiliary Valve

7



8MFA2006PTO






Y8MFA2005

Y8MFA2005PT2

WITH (2) 50 PSI AUXILIARY VALVES AND ADAPTER		
Part #	Adapter	Price
Y8MFA2001 PT2	None - 1/8" NPT Female	\$ 116.00
Y8MFA2002 PT2	Diss Female	119.00
Y8MFA2003 PT2	Diss Female w/hand nut	120.00
Y8MFA2006 PT2	Chemetron (NCG)	151.00
Y8MFA2005 PT2	Ohio (Ohmeda)	149.00
Y8MFA2008 PT2	Puritan	127.00
Y8MFA200M PT2	Medstar	170.00
Y8MFA2007 PT2	Oxequip	128.00

15 LPM DUAL AIR FLOWMETER

FIXED FLOW OXYGEN FLOWMETERS (1/2 - 15 LPM)



It's Durable - The lack of a flow tube and the durable metal body prevents breakage problems long associated with standard flowmeters.

It's Compact - The compact size allows the unit to fit into tight headwall areas and makes it virtually damage proof when other equipment is installed.

It's Accurate - The patented multi-spring technology insures accuracy at every flow setting without the need for a visible float as in Thorpe tube flowmeters

It's User-Friendly - The knurled flow adjustment knob and the large, easy-toread numbers make flow adjustment easy.

It's Versatile - The unit is available for headwall outlets as well as for ceiling mount drops. Both air and oxygen models are available.

FOR MORE FLOWMETERS, PLEASE VISIT OUR WEBSITE AT: WWW.RCMEDICAL.COM



	Adapters	Side Connection	Top Connection
	1/8" NPT FE	7MFA1001 \$68.00	7MFA1001TC \$68.00
	DISS Female Hand Tight	7MFA1003 \$74.00	7MFA1003TC \$74.00
	Ohmeda	7MFA1005 \$106.00	7MFA1005TC \$106.00
f	Chemetron (NCG)	7MFA1006 \$106.00	7MFA1006TC \$106.00
	Puritan-Bennett	7MFA1008 \$87.00	7MFA1008TC \$87.00

NEB BLOCK W/FLOWMETER

Features:

- The Neb Block allows a gas outlet to simultaneously or individually provide various respiration treatments without disconnecting any tubing or equipment.
- Provides easy to use On/Off switch to quickly access oxygen flow to a hand held nebulizer or resuscitation bag.
- Reliable fixed flow orifices provide six ranges of flow for aerosol therapy or resuscitation.

6-7 LPM 8-9 LPM 14-16 LPM 7-8 LPM 10-11 LPM

	6-7 LPM	8MFA6065	168.00
	7-8 LPM	8MFA6075	168.00
	8-9 LPM	8MFA6085	168.00
Ohmeda	10-11 LPM	8MFA6105	168.00
Offificad	14-16 LPM	8MFA6155	168.00
	6-7 LPM	8MFA6066	168.00
	7-8 LPM	8MFA6076	168.00
(D)	8-9 LPM	8MFA6086	168.00
12	10-11 LPM	8MFA6106	168.00
Chemetron (NCG)	14-16 LPM	8MFA6156	168.00
	6-7 LPM	8MFA6063	168.00
	7-8 LPM	8MFA6073	168.00
	8-9 LPM	8MFA6083	168.00
DISS Female Hand Tight	10-11 LPM	8MFA6103	168.00
	14-16 LPM	8MFA6153	168.00
	6-7 LPM	8MFA6068	168.00
	7-8 LPM	8MFA6078	168.00
VIII.	8-9 LPM	8MFA6088	168.00
Puritan Bennett	10-11 LPM	8MFA6108	168.00
	14-16 LPM	8MFA6158	168.00

SPECIALTY OXYGEN FLOWMETERS



15 LPM O2 Flowmeter for Chemetron Wall Outlet with 50 PSI-Chemetron Female Quick Connect

SFO-5 \$227.00



15 LPM O2 Flowmeter for Chemetron Wall Outlet with 50 PSI-Chemetron Female Quick Connect

SFO-2 \$249.00



Dual 15 LPM Oxygen Flowmeter with Chemetron Quick Connect Inlet and (2) Chemetron Couplers (Auxiliary Valves)

SFO-12 \$339.00



15 LPM O2 Flowmeter for Ohmeda Wall Outlet with 50 PSI-Ohmeda Female Quick Connect

SFO-4 \$227.00



15 LPM O2 Flowmeter for Ohmeda Wall Outlet with 50 PSI-Ohmeda Female Quick Connect

SFO-1 \$249.00



Dual 15 LPM Oxygen Flowmeter with Ohmeda Quick Connect Inlet and (2) Ohmeda Couplers (Auxiliary Valves)

SFO-11 \$339.00

SMALL CYLINDER CARTS (D&E)

#6105.....\$57.00



#PX-1060.....\$65.00



#6103-4.....\$222.00

41" high 3" casters swivel w/brake



These four-wheeled units offer better control and easy maneuverability in and out of tight spaces. They require no tilting or balancing, giving you total control over the cylinder load.

#1060-F \$177.00



Holds 4

#6205 **\$143.00**



Holds 2

LARGE CYLINDER CARTS AND STANDS (H,K)



Part #	Description	Price
6114	Single (Most Popular)	\$386.00
6114-PN	Single (Pneumatic Wheels)	418.00
6214	Dual	468.00
6214-PN	Dual (Pneumatic Wheels)	486.00

2 Wheels Only



Part # HCC2-M \$232.00



610FW-E.....\$276



Fiberglass Cylinder Bases 65045......\$98

SMALL CYLINDER CARTS (D&E)

2-Wheel Cart



	Cylinders		Dimensions
Part #	Held	Price	(Height-Depth-Width)
6041	4	261.00	42" - 10" - 18"
6061	6	242.00	42" - 15" - 18"
6121	12	318.00	42" - 20" - 22"

	With Sound Deadening Kit (below)					
	6041-SDK 4 312.00 42" - 10" - 18"					
	6061-SDK	6	296.00	42" - 15" - 18"		
6121-SDK 12 357.00 42" - 20" - 22"						

4-Wheel Cart



	Cylinders		Dimensions
Part #	Held	Price	(Height-Depth-Width)
6064	6	348.00	42" - 15" - 18"
6124	12	370.00	42" - 20" - 22"
6184*	18	435.00	42" - 30" - 16"
6244*	24	455.00	42" - 30" - 20"

With Sound Deadening Kit (below)						
6124-SDK 12 396.00 42" - 20" - 22"						
18	477.00	42" - 30" - 16"				
6244-SDK* 24 497.00 42" - 30" - 20"						
	12 18	12 396.00 18 477.00				

^{* -4&}quot; wheels

Multiple D & E Size Cylinder Stands With Wheels



	Cylinders		Dimensions
Part #	Held	Price	(Height-Depth-Width)
6040-PC	4	209.00	23" - 13" - 10"
6060-PC	6	231.00	23" - 12" - 15"
6120-PC	12	252.00	23" - 15" - 22"

With Sound Deadening Kit (below)					
6040-PC-SDK	4	244.00	23" - 13" - 11"		
6060-PC-SDK	6	267.00	23" - 12" - 15"		
6120-PC-SDK	12	282.00	23" - 15" - 22"		

ALL CARTS ABOVE ARE PAINTED GREEN

RED CARTS - ADD A "R" TO END OF PART NUMBER, PRICE IS \$65.00 MORE YELLOW CARTS - ADD A "Y" TO END OF PART NUMBER, PRICE IS \$65.00 MORE

SOUND DEADENING KITS (for Carts Above)



Anthony has developed an exclusive sound deadening system for their multiple cylinder carts, used to distribute cylinders in and around hospitals, rest homes, convalescent homes, etc. or where QUIET is essential.

These sound deadening kits include heavy-duty sound deadening matting and plastic bumper guards. This formula eliminates the rattling & clanging of steel & aluminum cylinders while in transport.

SOUND DEADENING KITS ALSO SOLD SEPARATELY PLEASE VISIT OUR WEBSITE: WWW.RCMEDICAL.COM





10" x 14" Plastic Sign Part # SP-51 \$37.00



10" x 14" Plastic Sign Part # SP-52 \$37.00

SMALL CYLINDER STANDS AND HOLDERS (D&E)

Multiple D & E Size Cylinder Stands



Part #	Cylinders Held	Price	Dimensions (Height-Depth-Width)
6040	4	184.00	20" - 13" - 11"
6060	6	197.00	20" - 12" - 15"
6120	12	209.00	20" - 15" - 22"
6180	18	251.00	20" - 15" - 32"
6240	24	251.00	20" - 20" - 32"
6400	40	403.00	20" - 40" - 25"

With	Sound	Deade	nina	Kit ((below)	
VVILII	Juliu	Deaue	IIIIIU	IXIL I	DCIUW	

6040-SDK	4	215.00	20" - 13" - 11"
6060-SDK	6	239.00	20" - 12" - 15"
6120-SDK	12	255.00	20" - 15" - 22"
6180-SDK	18	335.00	20" - 15" - 32"
6240-SDK	24	335.00	20" - 20" - 32"
6400-SDK	40	533.00	20" - 40" - 25"

Can be bolted to floor or free standing.

ALL STANDS ABOVE ARE PAINTED GREEN
RED STANDS - ADD A "R" TO END OF PART NUMBER, PRICE IS \$65.00 MORE
YELLOW STANDS - ADD A "Y" TO END OF PART NUMBER, PRICE IS \$65.00 MORE



MULTIPLE CYLINDER STANDS

CHROME



PX-1025

PX-1025 "D/E" Rack / 6 cylinders 9" x 15" x 10"\$ 134.00 PX-1027 "D/E" Rack / 12 cylinders 9" x 15" x 20" 190.00

Bolts to floor or free standing.



For any solid surfaces, crash carts, more...

Part # 610WM

\$97.00



For wheelchairs Part # 610WC \$99.00

FOR MORE INFORMATION ON SMALL CYLINDER BASES,

CYLINDER HOLDERS FOR BEDS, IV POLES, M6 CYLINDER CARTS AND STANDS,

AND MORE, PLEASE VISIT OUR WEBSITE AT:

WWW.RCMEDICAL.COM

OXYGEN REGULATORS - LITERS PER MINUTE

Responsive Respiratory Fixed Flow Oxygen Regulator CGA 870 for Small Cylinder (D & E)



Part #120-3040 \$47.00 for small cylinder (CGA 870)



Part #120-1061C \$47.00 for large cylinder (CGA 540)

ALL BRASS

0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 settings



Part #21017.....\$154.00 for small cylinder (870)

0, 1, 2, 4, 6, 8, 10, 12, 15, 18, 20, 25 settings



Two High Pressure **Outlets** Also

Part #21025.....\$169.00 1 - 25 LPM

Western Medica Oxygen Regulators

(Diaphragm Type)



Single stage regulator for large cylinders. 15 liters per minute. ORDER #M1-540-15 FG \$248.00

With Gauge-Savers ORDER #M1-540-15 FGP \$279.00

6 Year Warranty



Single stage flowmeter/regulator for large cylinders. 15 liters per minute.

ORDER #M1-540-15 FM....... \$276.00

Single stage flowmeter/regulator for small cylinders.

15 liters per minute.

ORDER #M1-870-15 FM...... \$285.00

MADE IN USA



Single stage regulator for small cylinders. 15 liters per minute. ORDER #M1-870-15 FG\$248.00

With Gauge-Savers ORDER #M1-870-15 FGP\$279.00



Single stage, swivel yoke regulator for small cylinders.

15 liters per minute.

ORDER #MR-870-15 FG \$246.00

With Gauge-Savers

ORDER #MR-870-15 FGP \$271.00

To order 8 or 5 liter regulators, add an 8 or a 5 to the end of our part number. The same price as our 15 liter regulators.

HIGH PRESSURE REGULATORS



Adjustable 0 - 100 psi with DISS Outlet

Gas	Part #	Cylinder Size	Price
Оху.	MI 540 PG	CGA 540 (Large)	\$258.00
Оху.	MI 870 PG	CGA 870 (Small)	\$241.00
Air	MI 346 PG	CGA 346 (Large)	\$294.00
Air	MI 950 PG	CGA 950 (Small)	\$274.00



Preset 50 psi with DISS Outlet

Gas	Part #	Cylinder Size	Price
Оху.	MIA 540 P	CGA 540 (Large)	\$241.00
Air	MIA 346 P	CGA 346 (Large)	\$261.00
N ₂ 0	MIA 326 P	CGA 326 (Large)	\$208.00
CO ₂	MIA 320 P	CGA 320 (Large)	\$235.00



Preset 50 psi with DISS Outlet

Gas	Part #	Cylinder Size	Price
Оху.	MIA 870 P	CGA 870 (Small D&E)	\$198.00
Air	MIA 950 P	CGA 950 (Small D&E)	\$259.00
N ₂ 0	MIA 910 P	CGA 910 (Small D&E)	\$208.00
CO ₂	MIA 940 P	CGA 940 (Small D&E)	\$254.00





Preset 50 psi with Ohmeda Quick Connect

Gas	Part #	Cylinder Size	Price
Оху.	SRO-19	CGA 540 (Large)	\$281.00
Air	SRO-34	CGA 346 (Large)	\$282.00



Preset 50 psi with Chemetron Quick Connect

Gas	Part #	Cylinder Size	Price
Оху.	SRO-21	CGA 540 (Large)	\$281.00
Air	SRO-36	CGA 346 (Large)	\$282.00



Preset 50 psi with Ohmeda Quick Connect

Gas	Part #	Cylinder Size	Price
Оху.	SRO-30	CGA 870 (Small D&E)	\$332.00
Air	SRO-32	CGA 950 (Small D&E)	\$284.00



Preset 50 psi with Chemetron Quick Connect

Gas	Part #	Cylinder Size	Price
Оху.	SRO-29	CGA 870 (Small D&E)	\$332.00
Air	SRO-31	CGA 950 (Small D&E)	\$284.00

SPECIALTY OXYGEN REGULATORS



Fixed flow Regulator With Ohmeda 50 psi Outlet. .5 -25 lpm 0 - 25 LPM Part # SRO-11..... \$305.00



Fixed flow Regulator With Chemetron 50 psi Outlet. .5 - 25 lpm Part # SRO-12..... \$305.00



Fixed flow Regulator With Puritan-Bennett 50 psi outlet. .5 - 25 lpm Part # SRO-13..... \$305.00



Oxygen Regulator/ Flowmeter 15 LPM, Large Cylinder with 50 PSI Valve Part # SRO-14......\$327.00



Oxygen Regulator with Dual Chemetron Female Couplers, 50 PSI, Large Cylinder Part # SRO-15..... \$337.00



Oxygen Regulator with Dual Ohmeda Female Couplers, 50 PSI, Large Cylinder Part # SRO-16..... \$337.00



Large Cylinder Oxygen Regulator with (2) 15 LPM Flowmeters Part # SRO-17...... \$331.00



Large Cylinder Oxygen Regulator with (2) 15 LPM Flowmeters and (2) Auxiliary Valves Part # SRO-18..... \$397.00



Oxygen Regulator, Preset 50 PSI with Ohmeda Quick Connect, Large Cylinder Part # SRO-19..... \$281.00



Oxygen Regulator, Preset 50 PSI with Puritan Bennett Quick Connect, Large Cylinder Part # SRO-20...... \$281.00



Oxygen Regulator, Preset 50 PSI with Chemetron Quick Connect, Large Cylinder Part # SRO-21...... \$281.00



Oxygen Regulator with 15 LPM flowmeter and 50 PSI Ohmeda Quick Connect, Large Cylinder Part # SRO-22..... \$364.00



Oxygen Regulator with 15 LPM Flowmeter and 50 PSI Chemetron Quick Connect, Large Cylinder Part # SRO-23...... \$364.00



Oxygen Regulator, Preset 50 PSI with Chemetron Quick Connect, Small Cylinder (CGA 870) Part # SRO-29..... \$332.00



0 - 25 LPM Regulator with Chemetron Female 50 PSI Outlet for Large Cylinder (CGA 540) Part # SRO-25..... \$254.00



0 - 25 LPM Regulator with Ohmeda Female 50 PSI Outlet for Large Cylinder (CGA 540) Part # SRO-26..... \$254.00

Specialty AIR Regulators, Visit our Website at: WWW.RCMEDICAL.COM

For Regulator Accessories See Page 31

HELIOX, AIR, CO₂ REGULATORS

Heliox Regulators



Part #	Ratio	For Cylinder size	Price
2210-20	70 - 30	Small D and E CGA 890 Yoke	\$297
2211-20	80 - 20	Small D and E CGA 890 Yoke	\$297
2212-20	70 - 30	Large CGA 280 Nut/Nipple	\$289
2213-20	80 - 20	Large CGA 280 Nut/Nipple	\$289

WESTERN MEDICAL REGULATORS (AIR or CARBON DIOXIDE)

Single-stage regulators preset to 50 psi or adjustable to 100 psi, or with a flowmeter. A sensitive neoprene diaphragm assembly provides optimum accuracy and sensitivity. Dual inlet filtering. Available for small or large cylinders.

	AIR	CO ₂
Preset	M1A-346P \$261.00	M1-320P \$235.00
50 PSI	for large cylinder	for large cylinder
Adjustable	M1-346PG \$294.00	M1-320PG \$379.00
0 - 100 PSI	for large cylinder	for large cylinder
12 Liters Per	M1-346-15FM \$480.00	M1-320-12FM \$379.00
Minute (lpm)	for large cylinder	for large cylinder
Preset	M1A-950P \$259.00	M1-940P \$254.00
50 PSI	for small cylinder	for small cylinder
Adjustable	M1-950PG \$274.00	M1-940PG \$330.00
0 - 100 PSI	for small cylinder	for small cylinder
12 Liters Per	M1-950-15FM \$499.00	M1-940-12FM \$390.00
Minute (lpm)	for small cylinder	for small cylinder

EASTERN RAIL SYSTEMS

Simple Outlet Expansion allows your hospital to use its existing outlets, and to create the number and placement of outlets needed. Simply mount a length of our rail on the wall - attach the needed outlets to the rail with our thumb screw clamp - attach a length of high-pressure hose to the block, and run it back to your existing outlet.





- Easy external mounting no need to shut down the unit.
- No more weight on outlets
- Shift the weight of humidifiers & suction bottles to the rail
- Eliminate overloading
- Reduce outlet leakage. Save time and money
- Increase efficiency
- Works with all piping systems
- Have outlets where you need them NOW!

Kit comes with:

1 each 2 Foot Rail
2 each Mounting Hardware
1 each 3 Foot Hose with Wall Adapter
1 each Double Gas Block

Outlet			
connection	GAS	Part #	Price
DISS	Oxygen	ERO-1DS	629.00
Chemetron	Oxygen	ERO-1CH	747.00
Ohmeda	Oxygen	ERO-10M	747.00
DISS	AIR	ERA-1DS	629.00
Chemetron	AIR	ERA-1CH	747.00
Ohmeda	AIR	ERA-10M	747.00
DISS	VAC	ERV-1DS	629.00
Chemetron	VAC	ERV-1CH	747.00
Ohmeda	VAC	ERV-1OM	747.00



FOR MORE EASTERN RAIL KITS

AND ACCESSORIES,

PLEASE GO TO OUR WEBSITE AT:

WWW.RCMEDICAL.COM



CUSTOM HOSE ASSEMBLIES

Ordering is easy as 1, 2, 3 . . . 4

How to Order:

Western Medica manufactures custom hose assemblies that conform to CGS, and NFPA standards. All assemblies are pressure-tested. Select your requirement from the charts below.

Example: OXY-024-FNN-OHF

1) OXY	Oxygen Service
2) 024	24" Long
3) FNN	Inlet Connection
4) OHF	Outlet Connection

Evacuation Hose (purple) price may be more, depending on your adaptor

3 Inlet Outlet Connection 1 Gas Service Connection 2 Hose Length

(Pricing for the above example)

FNN	22.00
OHF	80.00
Hose (.35 per inch x 24")	8.40
Assembly Charge	15.00
	\$125.40

MINIMUM HOSE LENGTH — 1 FOOT HOSE ORDERS FOR CUSTOM HOSES MAY NOT BE CANCELLED

Gas Service

OXY	Oxygen	Green
AIR	Medical Air	Yellow
SUC	Suction	White
N20	Nitrous Oxide	Blue
EVC	Evacuation	Purple
CO ₂	Carbon Dioxide	Gray
NIT	Nitrogen	Black
GNR	General	Clear

Adapters pictured below may be ordered Separately.

When ordering adapters only, please specify gas. Example: chm - oxy

All hose assemblies are 1/4" conductive (unless specified)



#CHM \$41.00 Chemetron (NCG) Male



#FNN \$22.00 Diss Female with Hex nut



#FHT \$24.00 Diss Female with Hand Nut



#OHM \$46.00 Ohmeda Male



#MBA \$24.00 Diss Male



#SHM \$26.00 Schrader Male



#OXM \$29.00 Oxequip Male



#PBM \$43.00 Puritan - Bennett Male



#CHF \$75.00 Chemetron (NCG) Female



#SHF \$114.00 Schrader Female



#PBF \$69.00 Puritan - Bennett Female



#OHF \$80.00 Ohmeda Female



#DBA \$25.00 Demand Body Adapter



#FP2 \$10.00 Female Pipe 1/8"



#MP2 \$10.00 Male Pipe 1/8"



#FP4 \$10.00 Female Pipe 1/4"



#MP4 \$10.00 Male Pipe 1/4"



#BLK \$0.00 Blank (No Connection)

REINFORCED MEDICAL HOSE



Hose only! \$5.00 per foot

Conductive, Non-UL Listed 400 psi Working Pressure 1500 psi Burst Pressure

MH4-04C	1/4" ID, Nitrous Oxide (Blue)
MH4-08C	1/4" ID, Carbon Dioxide (Grey)
MH4-12C	1/4" ID, Nitrogen (Black)
MH4-16C	1/4" ID, Medical Air (Yellow)
MH4-22C	1/4" ID, Suction (White)
MH4-24C	1/4" ID, Oxygen (Green)
MH5-22C	5/16" ID Suction (White)

FERRULE





For 1/4" hose .525 ID #M-625-A.....**\$1.00**

VICE-GRIPS

For 3/16", 1/4" and 5/16" ID Hose # C-5A **\$70.00**



COMMON HOSE ASSEMBLIES

DISS FEMALE HANDTIGHT X DISS FEMALE HANDTIGHT



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-100	\$74.00
5 Foot Air	CH-101	\$74.00
10 Foot Oxygen	CH-102	\$95.00
10 Foot Air	CH-103	\$95.00
15 Foot Oxygen	CH-104	\$116.00
15 Foot Air	CH-105	\$116.00
20 Foot Oxygen	CH-106	\$137.00
20 Foot Air	CH-107	\$137.00

DISS FEMALE HANDTIGHT X CHEMETRON MALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-108	\$91.00
5 Foot Air	CH-109	\$91.00
10 Foot Oxygen	CH-110	\$112.00
10 Foot Air	CH-111	\$112.00
15 Foot Oxygen	CH-112	\$133.00
15 Foot Air	CH-113	\$133.00
20 Foot Oxygen	CH-114	\$154.00
20 Foot Air	CH-115	\$154.00

DISS FEMALE HANDTIGHT X CHEMETRON FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-116	\$122.00
5 Foot Air	CH-117	\$122.00
10 Foot Oxygen	CH-118	\$143.00
10 Foot Air	CH-119	\$143.00
15 Foot Oxygen	CH-120	\$164.00
15 Foot Air	CH-121	\$164.00
20 Foot Oxygen	CH-122	\$185.00
20 Foot Air	CH-123	\$185.00

CHEMETRON MALE X CHEMETRON FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-124	\$140.00
5 Foot Air	CH-125	\$140.00
10 Foot Oxygen	CH-126	\$161.00
10 Foot Air	CH-127	\$161.00
15 Foot Oxygen	CH-128	\$182.00
15 Foot Air	CH-129	\$182.00
20 Foot Oxygen	CH-130	\$203.00
20 Foot Air	CH-131	\$203.00

DISS FEMALE HANDTIGHT X PURITAN BENNETT MALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-164	\$77.00
5 Foot Air	CH-165	\$77.00
10 Foot Oxygen	CH-166	\$98.00
10 Foot Air	CH-167	\$98.00
15 Foot Oxygen	CH-168	\$119.00
15 Foot Air	CH-169	\$119.00
20 Foot Oxygen	CH-170	\$140.00
20 Foot Air	CH-171	\$140.00

DISS FEMALE HANDTIGHT X OHMEDA MALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-132	\$92.00
5 Foot Air	CH-133	\$92.00
10 Foot Oxygen	CH-134	\$113.00
10 Foot Air	CH-135	\$113.00
15 Foot Oxygen	CH-136	\$134.00
15 Foot Air	CH-137	\$134.00
20 Foot Oxygen	CH-138	\$155.00
20 Foot Air	CH-139	\$155.00

DISS FEMALE HANDTIGHT X OHMEDA FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-140	.\$124.00
5 Foot Air	CH-141	.\$124.00
10 Foot Oxygen	CH-142	.\$145.00
10 Foot Air	CH-143	.\$145.00
15 Foot Oxygen	CH-144	.\$166.00
15 Foot Air	CH-145	.\$166.00
20 Foot Oxygen	CH-146	.\$187.00
20 Foot Air	CH-147	.\$187.00

OHMEDA MALE X OHMEDA FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-148	\$143.00
5 Foot Air	CH-149	\$143.00
10 Foot Oxygen	CH-150	\$164.00
10 Foot Air	CH-151	\$164.00
15 Foot Oxygen	CH-152	\$185.00
15 Foot Air	CH-153	\$185.00
20 Foot Oxygen	CH-154	\$206.00
20 Foot Air	CH-155	\$206.00

PURITAN BENNETT MALE X PURITAN BENNETT FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-156	\$120.00
5 Foot Air	CH-157	\$120.00
10 Foot Oxygen	CH-158	\$141.00
10 Foot Air	CH-159	\$141.00
15 Foot Oxygen	CH-160	\$162.00
15 Foot Air	CH-161	\$162.00
20 Foot Oxygen	CH-162	\$183.00
20 Foot Air	CH-163	\$183.00
	5-Foot Oxygen 5 Foot Air 10 Foot Oxygen 10 Foot Air 15 Foot Oxygen 15 Foot Air 20 Foot Oxygen	5-Foot Oxygen CH-156 5 Foot Air CH-157 10 Foot Oxygen CH-158 10 Foot Air CH-159 15 Foot Oxygen CH-160 15 Foot Air CH-161 20 Foot Oxygen CH-162

DISS FEMALE HANDTIGHT X PURITAN BENNETT FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-172	\$116.00
5 Foot Air	CH-173	\$116.00
10 Foot Oxygen	CH-174	\$137.00
10 Foot Air	CH-175	\$137.00
15 Foot Oxygen	CH-176	\$158.00
15 Foot Air	CH-177	\$158.00
20 Foot Oxygen	CH-178	\$179.00
20 Foot Air	CH-179	\$179.00

HUDSON RCI DISPOSABLE ADAPTERS AND CONNECTORS



No. 1667.....\$196.00 Ventilator Monitoring **Adaptor Circuit**

- Allows monitoring of ventilatory weaning parameters in intubated
- Includes oxygen nipple adaptor, aerosol tee, 6" CORR-A-FLEX tube, 15mm/22mm adaptor, m/m one-way valve and m/f one-way valve
- 50 per case



No. 1637.....\$195.00 Thermometer with Tee Adapter

- Combines No. 1635 Thermometer and No. 1077 Aerosol Tee
- Packaged 50 to a case



No. 1743.....\$291.00 Valved Neb Tee

- 22mm O.D. x 18mm I.D. x 22mm I.D. connections
- Tee may be left in-line when nebulizer is removed
- Port accepts most small volume nebulizer

More Adapters on our website

RCMedical.com

• 30 per case



No. 1648.....\$82.00 Aerosol Tee Connector

- 22mm O.D. x 15mm I.D./22mm O.D. x 22mm O.D.
- 50 per case



No. 1641.....\$87.00 Elbow Adaptor

- Standard 22mm I.D. connector on vertical axis
- 22mm O.D. connection on horizontal axis
- Packaged 50 to a case



No. 1642.....\$122.00

Pressure Line Adaptor

- 22mm I.D./O.D. connections
- Barbed elbow accepts 5 to 7mm I.D. tubing and swivels 360 degrees
- Packaged 50 to a case



No. 1623.....\$145.00 Patient Elbow Adaptor

- 15mm O.D. x 15mm I.D./22mm O.D.
- 50 per case



No. 1644.....\$134.00 One-Way Valve

- Standard 22 mm I.D./O.D. connections
- Valve offers low flow resistance and
- a capped monitoring port
- Packaged 50 to a case



No. 1638.....\$87.00 **Aerosol Tee Connector**

- Horizontal ports of the tee are 22mm O.D.
- Vertical port is 18mm I.D./22mm O.D.
- Packaged 50 to a case



No. 1077.....\$64.00 Aerosol Tee Connector

- Anti-spill design minimizes aspiration of liquids
- Horizontal axis has standard 22mm O.D. connections
- Vertical axis incorporates a standard 15mm I.D. connection and two 22mm O.D. concentric rings
- Packaged 50 to a case



No. 1639.....\$87.00 Nebulizer Tee Connector

- Standard 22 mm I.D./O.D.
- connections
- Nebulizer port is 18mm I.D. which accepts most small volume nebulizers
- Packaged 50 to a case



No. 1421.....\$92.00 Universal Cuff Adaptor

- Standard 22mm I.D. connections
- Packaged 50 to a case



No. 1422.....\$98.00 Multi Adaptor

- Standard 15mm I.D./22mm O.D. connections
- Packaged 50 to a case



No. 1409.....\$69.00 Large Tubing Connector

- Standard 22mm O.D. connections
- Packaged 50 to a case



No. 1640.....\$87.00 **Tubing Adaptor**

- Standard 22mm I.D./O.D.
- Packaged 50 to a case



No. 1751.....\$290.00 MDI Adapter

- 15mm O.D x 15mm I.D. connections for in-line placement in 15mm tubing
- 50 per case



No. 1659.....\$244.00 MDI Adapter

- 22mm O.D x 22mm I.D. connections for in-line placement in 22mm tubing
- 50 per case



No. 1423.....\$131.00 **Tubing Adaptor**

- 22mm I.D. x 5 7mm tapered nipple
- Packaged 50 to a case



No. 1669.....\$196.00 Monitoring Adaptor

• 15mm O.D. x 15mm I.D./22mm

- O.D.
- Standard luer fitting
- 50 per case

DISPOSABLES



#2200 Flexible Connector (Silicone Rubber)

22mm I.D. x 22mm I.D.

This adapter was designed almost identical to the ones listed below:

INTEC 2200 MARQUEST 172200
BENNETT 3443 AM SYSTEMS 1139
Instrumentation Ind. 101-00-3443

Part Number 2200 \$294.00 per case of 50 \$9.00 each



Tapered Silicone Rubber Connector

27mm ID x 22mm ID

Part Number 1426 \$134.00 per case of 10



MDI Accessory by Hudson

The POCKET CHAMBER is equipped with a highly sensitive silicone valve that opens easily without restricting the air flow. This one-way valve helps keep patients from accidentally blowing medication out of the chamber, and may aid patients experiencing difficulty timing their inhalation.

Part Number 1001-10 \$196.00 per case of 10



Main Flow Bacterial / Viral Filter

- 22mm ID x 15mm ID/22mm OD
- Bacterial filtration efficiency 99.999+%
- Viral filtration efficiency 99.99+%
 Dead space: 42mL
 - Flow resistance: approximately
- 1.5 cmH₂O at 60 LPM
 Individually packaged 50 to a case

Part Number 1605 \$198.00 per case of 50



Supply Line Adapter

For extending smooth bore oxygen tubing Tapered for excellent gripping

> Part Number 65100 \$163.00 per case of 200



Corrugated Tubing

- Continuous 100-foot roll of tubing with cuttable 6-inch sections
- Packaged in protective, convenient, plastic-lined dispenser box

Part Number 1680 \$52.00 per roll



Ventilator Circuits

Adult Standard Ventilator Circuit, 72" long. Use on Puritan-Bennett 7200A; Servo 300, 900A, 900B, and 900C; Hamilton Veolar; Emerson; Engstrom Erica; PPG Irisa; Infrasonics Adult Star and Bird 8400ST

> Part Number 1613 \$149.00 per case of 10



Voldyne Incentive Spirometer

Voldyne volumetric exerciser features Good - Better - Best flow window in a compact, ergonomic design with advanced, low work-of-breathing filter design and built in handle. 5000 ml

> Part Number 8884719010 \$158.00 per case of 12



E-Z Wrap

Soft closed-cell foam tubes designed to prevent ear chafing. Light gray color Packed 50 pair per case.

> Part Number 1016 \$73.00 per case of 50



Noseclip

Part Number 64019 \$87.00 per case of 100



Aerosol Drainage Bags

Disposable Aerosol Drainage Bags with tee adaptor fits 22mm ID corrugated tubing. Bag has plastic chain for securing and convenient drainage port. 750 cc capacity. With wye adaptor.

> Part Number 1742 \$182.00 per case of 50

RESPIROMETERS

MARK 8

MARK 14

HALOSCALE STANDARD





tidal volumes. Two dials. Small dial records 0-1

liter. Large dial 0-100 liters.

1.4" Face

For measurement of minute volumes as well as Specially suited for the measurement of tidal volumes and may also be used for measurement of minute volumes. Two dials. Small dial records 0-10 liters. Large dial records 0-1 liter.



The original design of Haloscale Respirometer incorporating concentric dials and an enlarged face for improved readability. Equally convenient for measuring both tidal and minute volumes. Large 57 mm dial diameter. Red/black hands and black scales on a white background with yellow highlight ring. Two concentric outer dials. First dial records 0-1 liter. Second dial records 0-100 liters.

RESPIROMETERS

ORDER #	DESCRIPTION	PRICE
2750	Mark 8 Respirometer	\$3405.00
2751	Mark 14 Respirometer	3405.00
2768	Haloscale Standard	3455.00

RESPIROMETERS WITH SPIRO GUARD

2756	Mark 8 Respirometer	3469.00
2757	Mark 14 Respirometer	3469.00
	Haloscale Standard	

RESPIROMETERS WITH SPIRO GUARD AND EXTRA LARGE CARRYING CASE

2762	Mark 8 Respirometer	3499.00
2763	Mark 14 Respirometer	3499.00
2772	Haloscale Standard	3549.00

Used Respirometers

Fully Serviced, Cleaned, Oiled and Calibrated

1 Year Warranty

Order #2799 \$1675 (Case included)

ELECTRONIC DIGITAL RESPIROMETER

A respirometer long overdue in this modern electronic age

Part #



ANESTHESIA ASSOCIATES

295	Respirometer	\$1608.00
Options		
296	6' Retractable Cable (extended length)	\$185.00
298	Carrying Case	118.00
299	Spare Flow Sensing Unit	594.00

Model 295, a two piece detachable respirometer comes complete ready to use with 9v battery and one (1) 22 mm female x 22 mm female adapter.



Rubber Protective Cover for Small Faced Respirometers (mk. 8, 14)

Part # 2700 Spiro-Guard \$98.00

MODEL 295 DIGITAL RESPIROMETER IS AN ELECTRONIC VERSION OF THE MECHANICAL RESPIROMETER. IT DOES THE SAME JOB PLUS

Easier to read

· Saving on the initial purchase and maintenance

• More reliable

· Extremely durable

• No gears, shafts, or bearings to break down in the electronic unit

• Separates for ease of patient use, sterilization and cleaning

REPAIRS

Price

\$369.00

Includes cleaning, calibration and oiling. Parts are extra

All Models Repaired including: Wright • Ferraris • Haloscale **BOC • Fraser - Harlake • Ohmeda Boehringer**

NPT PIPE FITTINGS

Female Tee



	Part No.
1/8 "	NPT 7034
1/4"	NPT 7035
	440.00

\$10.00

Street Tee



	Part No.	
¹/8″ NPT	7032	
¹ / ₄ " NPT	7033	
\$10.00		

Cross



Part No.		
¹ / ₈ " NPT 7036		
¹ / ₄ " NPT 7037		
\$14.00		

45° Street Elbow



	Part No.
¹ / ₈ " NPT	7056
¹ / ₄ " NPT	7057
\$10.00	

Connectors



			Part No
1/8" Fe	emale x 1/8"	Male	1017
¹/8 ″ F€	emale x 1/4"	Male	7024
¹/₄″ F€	emale x 1/8"	Male	702
¹/₄″ F€	emale x 1/4"	Male	7026
\$10.00			

Hex Nipple



		Par	t No.
1/8 "	Male x 1/8"	Male	7020
1/4 "	Male x 1/4"	Male	7021
1/8 "	Male x 1/4"	Male	7022

\$10.00

Hose Barb Connector



Part No.				
1/8" Male x 1/8" [Dia 7038			
1/8" Male x 3/16" [Dia 7039			
1/8" Male x 1/4" [Dia 7040			
¹ / ₄ " Male x ¹ / ₄ " Dia 7043				
\$10.00				

Coupler



			Part No
1/8 "	Female x 1/8	Female	7027
1/4 "	Female x 1/4	Female	7028
1/4 "	Female x 1/8	Female	7029

\$10.00

90° Elbow



Part No. 1/8" Female x 1/8" Female **7058** 1/4" Female x 1/4" Female **7059**

\$10.00

Plug



	Part No
¹ /4" NPT	701
¹ /8" NPT	701
\$10	0.00

Hose Barb Connector



	Part No.	
¹ / ₈ " Female x ¹ / ₄ " Dia	7046	
¹ / ₄ " Female x ¹ / ₄ " Dia	7049	
\$10.00		

Cap



	Part No.	
¹ / ₄ " NPT	7016	
¹ / ₈ " NPT	7018	
\$10.00		

Pipe Tee



Part No. ¹/₈" NPT......**7030** ¹/₄" NPT......**7031** \$10.00

Drop Ear



Part No. Tee 1/8"F x 1/8"F x 1/8"F7081 Elbow 1/8" F x 1/8" F...... 7082 \$21.00

Pipe Nipple



	Part No.		
1/8" Male x 11/2" Long	7050		
1/8" Male x 2" Long	7051		
1/8" Male x 21/2" Long	7052		
1/4" Male x 11/2" Long			
1/4" Male x 2" Long			
1/4" Male x 21/2" Long			
1/8" Male x 31/2" Long			
\$14.00			

90° Elbow



	Part No.
¹ / ₈ "F x ¹ / ₈ "M	7063
¹ / ₄ "F x ¹ / ₄ "M	7065
\$10.00	

"Y" BLOCK ASSEMBLIES

Ohmeda Style

Chemetron Style

(NCG)

Ohmeda Style "Y" Assembly-"Stubby"

(Ohio)



2152	O ₂	193.00
2252	Air	193.00
2352	Vac	193.00



3152	O ₂	195.00
3252	Air	195.00
3352	Vac	195.00

Puritan-Bennett

Style

4152	O ₂	193.00
4252	Air	193.00
4352	Vac	193.00



2162	O ₂	201.00
2262	Air	201.00
2362	Vac	201.00

T BLOCKS



Gas Part # Price Oxy 2152-T \$203.00 Air 2252-T \$203.00 Vac 2352-T \$203.00



OHMEDA Style Gas Part # Price Oxy 2152-TL \$203.00 Air 2252-TL \$203.00 Vac 2352-TL \$203.00



 Gas
 Part #
 Price

 Oxy
 2152-TR
 \$203.00

 Air
 2252-TR
 \$203.00

 Vac
 2352-TR
 \$203.00

CHEMETRON Style



Part #	Price
3152-T	\$203.00
3252-T	\$203.00
3352-T	\$203.00
	3152-T 3252-T



	Price
152-TL	\$203.00
252-TL	\$203.00
352-TL	\$203.00
	Part # 152-TL 252-TL 352-TL



 Gas
 Part #
 Price

 Oxy
 3152-TR
 \$203.00

 Air
 3252-TR
 \$203.00

 Vac
 3352-TR
 \$203.00

TRI-BLOCK ASSEMBLIES

2

Ohmeda Style

SP-112	O ₂	295.00
SP-113	Air	295.00
SP-114	Vac	295.00



Chemetron Style

SP-115	O ₂	295.00
SP-116	Air	295.00
SP-117	Vac	295.00

DISS "Y" ASSEMBLY WITH CHECK VALVES



7211	O ₂	95.00
7212	Air	95.00
7213	Vac	95.00

FOR MORE "Y" BLOCK ASSEMBLIES
PLEASE VISIT OUR WEBSITE
WWW.RCMEDICAL.COM

DISS "Y" WITH 2 OHMEDA COUPLERS

DISS "Y" WITH 2 CHEMETRON COUPLERS



SP-214	O ₂	223.00
SP-214A	Air	223.00
SP-214V	Vac	223.00



SP-215	O ₂	223.00
SP-215A	Air	223.00
SP-215V	Vac	223.00

QUICK CONNECTS (MALE)			
	Ohmeda		Chemetron
1/8" NPT Male	Gas Part # Price Oxygen .2101 \$41.00 Air .2201 .41.00 Vacuum .2301 .41.00	1/8" NPT Male	Gas Part # Price Oxygen 3101 \$37.00 Air 3201 37.00 Vacuum 3301 37.00
1/8" NPT Female	Oxygen2102\$43.00 Air220243.00 Vacuum230243.00	1/8" NPT Female	Oxygen3102\$39.00 Air320239.00 Vacuum330239.00
1/4" NPT Male	Oxygen	1/4" NPT Male	Oxygen
DISS Male	Oxygen	DISS Male	Oxygen3105\$47.00 Air320547.00 Vacuum330547.00
DISS Hex Nut	Oxygen2106\$47.00 Air220647.00 Vacuum230647.00	DISS Hex Nut	Oxygen3106\$44.00 Air320644.00 Vacuum330644.00
DISS Hand Tight	Oxygen	DISS Hand Tight	Oxygen
1/4" Hose Barb	Oxygen	1/4" Hose Barb	Oxygen3111\$41.00 Air321141.00 Vacuum331141.00
DISS Male with Check Val	with Check Valve DISS Male with Check Valve		lve
	Oxygen 2114 \$50.00 Air 2214 50.00 Vacuum 2314 50.00		Oxygen
	Puritan	Bennett	
1/8" NPT Male	Gas Part # Price Oxygen .4101 \$37.00 Air .4201 .37.00 Vacuum .4301 .37.00	DISS Hex Nut	Gas Part # Price Oxygen 4106 \$41.00 Air 4206 41.00 Vacuum 4306 41.00
1/8" NPT Female		DISS Hand Tight	
	Oxygen .4102 \$39.00 Air .4202 .39.00 Vacuum .4302 .39.00		Oxygen
1/4" NPT Male	1/4" NPT Male 1/4" Hose Barb w/Knob		
	Oxygen .4103 \$41.00 Air .4203 .41.00 Vacuum .4303 .41.00		Oxygen4111\$43.00 Air421143.00 Vacuum431143.00
DISS Male with Knob		DISS Male w/Knob & C	
	Oxygen4105\$44.00 Air 4205 44.00		Oxygen

24 RC Medical

Air.....4214.....48.00

Vacuum 4314......48.00

Air.....4205......44.00

Vacuum 430544.00

QUICK CONNECTS (FEMALE)

CHEMETRON STYLE CHECK UNITS (Female Couplers)

DISS Male

DISS Female



Gas	Part No.
O ₂	3150
Air	3250
Vac*	3350



Gas	Part No.
O ₂	3155
Air	3255
Vac*	3355

Gas	Part No.
O ₂	3155
Air	3255
Vac*	3355

1/8" NPT Female

Gas	Part No.
O ₂	2150
Air	2250
Vac*	2350



OHMEDA STYLE CHECK UNITS (Female Couplers)

DISS Male

Gas	Part No.
O ₂	2155
Air	2255
Vac*	2355



Part No.
3153
3253
3353



Gas	Part No.
Hand Tight O2	3158
Hand Tight Air	3258
Hand Tight Vac*	3358

1/4" NPT Male



Gas	Part No.
O ₂	2153
Air	2253
Vac*	2353



Gas	Part No.
Hand Tight O2	2158
Hand Tight Air	2258
Hand Tight Vac*	2358

\$87.00

PURITAN-BENNETT STYLE CHECK UNITS (Female Couplers)

\$81.00



Gas	Part No.
O ₂	4150
Air	4250
Vac*	4350

DISS Male

DISS Female



Gas	Part No.
O ₂	4155
Air	
Vac*	4355

SCHRADER STYLE CHECK UNITS (Female Couplers)



Gas	Part No.
O ₂	6150
Air	6250
Vac*	6350

DISS Male

DISS Female



Gas	Part No.
O ₂	6155
Air	6255
Vac*	6355

1/4" NPT Male

1/8" NPT Female



Gas	Part No.
O ₂	4153
Air	4253
Vac*	4353

Gas	Part No.
Hand Tight O2	4158
Hand Tight Air	4258
Hand Tight Vac*	4358

1/4" NPT Male

1/8" NPT Female



Gas	Part No.
O ₂	6153
Air	6253
Vac*	6353

Gas	Part No.
Hand Tight O2	6158
Hand Tight Air	6258
Hand Tight Vac*	

\$117.00

CUSTOM QUICK CONNECTS AND COUPLERS

OHMEDA FEMALE X OHMEDA FEMALE



Gas	Part #	Price
Oxygen	SP-105	\$174.00
Air	SP-105A	174.00
Vacuum	SP-105V	174.00

OHMEDA FEMALE COUPLER X DISS FEMALE HAND TIGHT WITH 90° ELBOW



Gas	Part #	Price
Oxygen	SP-51-OHF	\$139.00
Air	SP-51-AOHF	139.00
Vacuum	SP-51-VOHE	139 00

OHMEDA OFFSET "V" BLOCK



Gas	Part #	Price
Oxygen	SP-211	\$237.00
Air	SP-211A	237.00
Vacuum	SP-211\/	237 00

CHEMETRON FEMALE X CHEMETRON FEMALE



Gas	Part #	Price
Oxygen	SP-106	\$188.00
Air	SP-106A	188.00
Vacuum	SP-106V	188.00

CHEMETRON FEMALE COUPLER X DISS FEMALE HAND TIGHT WITH 90° ELBOW



Gas	Part #	Price
Oxygen	SP-53CHF	\$139.00
Air	SP-53ACHF	139.00
Vacuum	SP-53V/CHF	139 00

CHEMETRON OFFSET "V" BLOCK



Gas	Part #	Price
Oxygen	SP-210	\$237.00
Air	SP-210A	237.00
Vacuum	SP-210V	237.00

CHEMETRON FEMALE COUPLER X OHMEDA MALE QUICK CONNECT



Gas	Part #	Price
Oxygen	SP-129	\$149.00
Air	SP-129A	149.00
Vacuum	SP-129V	149.00

CHEMETRON MALE QUICK CONNECT X OHMEDA FEMALE COUPLER



Gas	Part #	Price
Oxygen	SP-130	\$149.00
Air	SP-130A	149.00
Vacuum	SP-130V	149.00

PURITAN BENNETT MALE X OHMEDA FEMALE COUPLER



Gas	Part #	Price
Oxygen	SP-209	\$129.00
Air	SP-209A	129.00
Vacuum	SP-209V	129.00

PURITAN BENNETT MALE X CHEMETRON FEMALE COUPLER



Gas	Part #	Price
Oxygen	SP-208	\$129.00
	SP-208A	
Vacuum	SP-208V	129.00

OHMEDA VERTICAL "T" BLOCK



Gas	Part #	Price
Oxygen	SP-199	\$237.00
	SP-199A	
Vacuum	SP-199V	237.00

CHEMETRON VERTICAL "T" BLOCK



Gas	Part #	Price
Oxygen	SP-196	\$237.00
Air	SP-196A	237.00
Vacuum	SP-196V	237.00

RC MEDICAL CAN ASSEMBLE MOST CUSTOM CONFIGURATIONS
PLEASE EMAIL: QUESTIONS@RCMEDICAL.COM WITH DRAWINGS OR REQUIREMENTS

DISS FITTINGS



DISS 90° Elbow Male x Female

Oxygen	Air	Vac
#0705	#SP-93	#SP-94
\$23.00	\$42.00	\$42.00



DISS Male x NPT Male

Oxygen	Air	Vac
1/8"#0701\$11.00	#0731\$15.00	#0751\$15.00
1/4"#0703\$11.00	#0733\$15.00	#0753\$15.00



DISS Male x DISS Male

Oxygen	Air	Vac
#0704	#SP-91	#SP-92
\$12.00	\$42.00	\$42.00

DISS Male x 1/8" NPT Male with Check Valve



Oxygen	Air	Vac
1/8" #0712\$15.00	#0742\$14.00	#0762\$22.00
1/4" #0713\$15.00	#0743\$14.00	#0763\$22.00

DISS Female x DISS Female



Oxygen	Air	Vac
#0104	#SP-88	#SP-89
\$21.00	\$29.00	\$29.00

DISS Male x DISS Female



Oxygen	Air	Vac
#SP-95	#SP-96	#SP-97
\$35.00	\$35.00	\$35.00



DISS Female x 1/8" NPT Male

C	Oxygen	Air	Vac
	#0110	#0209	#0309
_	\$14.00	\$13.00	\$15.00

DISS Male x 1/8" NPT Female



Oxygen	Air	Vac
#0702	#SP-98	#SP-99
\$14.00	\$16.00	\$16.00



50 PSI Auxiliary Valve

OPTO - \$21.00 (OXYGEN)

APTO - \$21.00 (AIR)

Thread Comparisons Guide:





(Measurements given are approximate.)

For Oxequip, Medstar and Schrader quick connects

PLEASE VISIT OUR WEBSITE AT: WWW.RCMEDICAL.COM





NEB DRIVE RESCUE FLOW

An Economical adapter that delivers a preset flow from standard wall outlets. It's a convenient, dependable way to administer aerosols.

DISS FEMALE



\$36.00

O,	6 to 7 LPM	8116
O ₂	7 to 8 LPM	8117
O,	8 to 9 LPM	8118
O,	10 to 11 LPM	8119
Air	6 to 7 LPM	8216
Air	7 to 8 LPM	8217
Air	8 to 9 LPM	8218
Air	10 to 11 LPM	8219
Rescue Flow	14 to 16 LPM	8115
O ₂	24 to 26 LPM	8114

OHMEDA



\$55.00

O,	6 to 7 LPM	8126
O ₂	7 to 8 LPM	8127
O ₂	8 to 9 LPM	8128
O ₂	10 to 11 LPM	8129
Air	6 to 7 LPM	8226
Air	7 to 8 LPM	8227
Air	8 to 9 LPM	8228
Air	10 to 11 LPM	8229
Rescue Flow	14 to 16 LPM	8125
O ₂	24 to 26 LPM	8124

CHEMETRON



\$55.00

6 to 7 LPM	8136
7 to 8 LPM	8137
8 to 9 LPM	8138
10 to 11 LPM	8139
6 to 7 LPM	8236
7 to 8 LPM	8237
8 to 9 LPM	8238
10 to 11 LPM	8239
14 to 16 LPM	8135
24 to 26 LPM	8134
	7 to 8 LPM 8 to 9 LPM 10 to 11 LPM 6 to 7 LPM 7 to 8 LPM 8 to 9 LPM 10 to 11 LPM

PURITAN-BENNETT



\$41.00

O,	6 to 7 LPM	8146
O ₂	7 to 8 LPM	8147
O ₂	8 to 9 LPM	8148
O,	10 to 11 LPM	8149
Air	6 to 7 LPM	8246
Air	7 to 8 LPM	8247
Air	8 to 9 LPM	8248
Air	10 to 11 LPM	8249
Rescue Flow	14 to 16 LPM	8145
O,	24 to 26 LPM	8144

To order Neb Drive or Rescue Flow with shut off valve, please add an "S" to the end of the part number.

The price is \$40.00 more.



Shut off valve allows you to leave neb drive in wall outlet.

REPLACEMENT RESPIRATORY OXYGEN SENSORS

Replaces: Ohmeda 6050-0004-110



Part # SS-20A



Part # SS-03

OEM Compatible with: Analytical Ind.: PSR-11-77 C&S: G0-070 Ceramatec: CAG-17 City Technology: MOX-4 Datex-Ohmeda: 6051-0000-219 Envitec: OOM105 Hamilton: HM-03 Hudson RCI: V-03, 5558 Integral Process: 65013 Maxtec: MAX-17, R116P10 Mercury Medical: 10-103-02 MSA: 804674 Newport: OM 100, VM 200 Nuova: E-77/O Sensidyne: SV-03 SensorTech: ST-03 Teledyne: T-7 Ventrex: V-03 Vlamed: T-7 VTI: 103102



Part # SS-11A

OEM Compatible with:

Air Shields: 6735142
Alphamed: R-17
Analytical Ind.: PSR-11-917-J
Smiths Medical > BCI: R-17
BioMed Devices: M-2, M-10
Bird Products: 66030
C82: G0-030
Ceramates: CAG-13
Datex-Ohmeda: 6051-0000-217
Draeger: 6735142
Envitec: OOM103-1, OOM202-1
Hamilton: HM-11
Hudson RCI: 5567, V-11
Invivo: 9445
Lifecare: PLV-100, 37-079
Maxtec: MAX-13, R115P10
Mercury Medical: 10-103-05
MSA: 406931
Oxitron: 3310
Pacifitech: PT-11A



Part # SS-25A

OEM Compatible with: Analytical Ind.: PSR-11-917-J1 C&S: G0-250 Maxtec: MAX-250ESF, R115P85 Mercury Medical: 10-103-08 Nuova: E-25/0 Pacifitech: PT-25A Sensidyne: SV-25A VTI: 103104



Part # SS-12A

OEM Compatible with: Analytical Ind.: PSR-11-917M C&S: G0-040 Ceramatec: CAG-12 City Technology: MOX-9 Datex-Ohmeda: 6051-0000-222 Envitec: OOM202 Hamilton Medical: 396008, 396009, HM-12 Hudson RCI: 5566, V-12 Integral Process: 65016 IMT Medical: PE300 Maxtec: MAX-12, R109P45 Mercury Medical: 10-103-06 MSA: 472062, 470644 Nuova: E-17/M Pacifitech: PT-12A Medtronic > Puritan Bennett: 4-072214-00 Siemens: 9004979-E347E, 6419332-E380E Sensidyne: SV-12A SensorTech: ST-01 Spacelabs: 90513 Taema: YR049700 Teledyne: R-24 MED, B-62321 Ventrex: V-12A Vlamed: R-24Med VTI Vascular Technology: 103100

ALL SENSORS \$149.00

Respironics: 37-079 Sensidyne: SV-11A

SensorTech: ST-11

Teledyne: R-17 Vlamed: R-41V VTI: 100580

Unconditional 12 month Warranty!
Exceptional Quality and Performance!
Canned for increased shelf life and freshness!



FOR MORE SENSORS
PLEASE VISIT OUR WEBSITE
WWW.RCMEDICAL.COM

MINIOX 200E OXYGEN ANALYZER



MiniOX 200E Oxygen Analyzer

Part # 600210......\$914.00
Comes complete with sensor, cable, tee and AA batteries
Replacement Sensor # 606940.....\$179.00

The MiniOX 200E is best used anywhere an inline analyzer is needed, when alarms are not deemed important. It is especially useful in conjunction with older ventilators or oscillators. Designed with simplicity and portablility in mind, this analyzer quickly adapts for use in tubing applications and is ideal for checking oxygen percent accuracy in line or on concentrations. Its patented ergonomic design fits comfortably in your hand and features one-touch on/off and calibration. The MiniOX 200E will calibrate to room air (20.9%) or 100% Oxygen. It is ideal for applications needing an external sensor.

1

2

- Large LCD Easy-to-read characters
- Calibration Reminder Reminds once a week
- 3 One-Touch Calibration Auto calibrate to room air or 100% Oxygen
- (4) On/Off Button One-touch on and off
- 5 Patented Ergonomic Design Fits comfortably in your hand
- 6 External Sensor Cable 10 ft. extended, 3 ft. retracted
- 7 Low Battery Indicator 5,000 + Hours of battery life
- (8) Rubber Gasket Prevents moisture from contact board
- External Sensor MiniOX 250E, 0-100% Oxygen, 24 month warranty

Approximately 5,000 hours with continuous use

6

MISCELLANEOUS

PLASTIC BAGS (DISPOSABLE)





PATIENT SET-UP (Cotton Drawstring)

Part No.	Description	Quantity	Price
0049	8 x 10	500	\$131.00
0059	11 3/4 x 15	500	\$157.00
XL-S2	16 x 18 x 3	250	\$148.00

EQUIPMENT COVERS

Dust covers designed to protect ventilators and other equipment. Comes on a roll.

Part No.	Description	Quantity	Price
PLC-60	28 x 22 x 58 (clear)	55	\$157.00
PLC-60/B	28 x 22 x 58 (buff)	55	\$157.00
0084	28 x 22 x 56 (clear***)	50	\$127.00
7083	28 x 22 x 56 (yellow)	50	\$115.00
0126	48 x 25 x 42 (red)	50	\$147.00

^{***}printed, "CLEAN, READY FOR USE"

Main Flow Bacterial / Viral Filter




- 22mm ID x 15mm ID/22mm OD
- Bacterial filtration efficiency 99.999+%
 - Viral filtration efficiency 99.99+%
 - Dead space: 42mL
 - Flow resistance: approximately
 1.5 cmH₂O at 60 LPM
 - Individually packaged 50 to a case

Part Number 1605 \$198.00 per case of 50

ANESTHESIA ASSOCIATES PRESSURE MANOMETER



Part#	Description	Price
199	Gauge with 1/4" NPT	
274	Gauge with 5-55 cm H ₂ C adjustable pressure relief valve and "T" for supply tubing.	244.00



Manometer has range of -40 to +80 cm H.O in a 2½" diameter chrome plated case with ½" NPT bottom connection. Has front zero correcting screw.



MAX is a resetable low friction pointer to capture maximum pressure or vacuum reading

POSEY CUFFLATOR™

ENDOTRACHEAL TUBE INFLATOR & MANOMETER



An accurate, convenient tracheal cuff inflator and manometer. No more syringes, electricity or complicated connections needed to check cuff pressure regularly. Air vent button and inflator bulb let you vary cuff pressure quickly and easily. Clip on back allows easy positioning. Gauge shows recommended pressure range in cm H₂O.

#8199 Cufflator with instrument case......\$679.00

FOR MORE PRODUCTS AND SPECIALS PLEASE VISIT OUR WEBSITE

WWW.RCMEDICAL.COM

#2200 Flexible Connector (Silicone Rubber)

22mm I.D. x 22mm I.D.

This adapter was designed almost identical to the ones listed below:

INTEC 2200 MARQUEST 172200
BENNETT 3443 AM SYSTEMS 1139
Instrumentation Ind. 101-00-3443



Part Number 2200 \$294.00 per case of 50 \$9.00 each

29

MISCELLANEOUS

Timeter Airdyne Hospital Grade Air Compressors



Uniquely designed for bed-side use. Special material absorbs most noise for quiet performance. Audible and visual high temperature and low pressure alarms with DC powered six-volt backup. Vibration-free gauge shows diminished pressure, impaired drying or exceeded flow capacity. Circuit breaker and thermal overload protector included. Continuous flow of 64 LPM at 50 psi.

Part Number TE20 \$4823.00



Designed for bedside use. Quiet, vibration-free performance. Provides condensation-free compressed air to respiratory therapy equipment. Vibration-free gauge shows diminished pressure, impaired drying or exceeded flow capacity. Circuit breaker and thermal overload protector included. Continuous flow of 42 LPM at 50 psi.

Part Number TE35 \$4113.00

DeVilbiss Pulmo-Aide® Medication Nebulizer



3 PRONG HOSPITAL GRADE PLUG

Durable AC-powered medication nebulizer with Pulmo-Neb® 60-day disposable nebulizer and tubing. Unit weighs 7 lbs. and has 9 - 16 LPM range. 5 year warranty. **115 vac**

DV4650D.....\$180.00

BLENDERS - PRECISION MEDICAL

Air/O



Its unique. It is the only unit available that allows the clinician to turn the bleed control to ON or OFF providing quieter operation. This feature is especially useful in the NICU and PICU.

Part #	Flow	Rate	Price
5200	High Flow	15 - 120 lpm	999.00
5300	Low Flow	3 - 30 lpm	999.00

Heliox



Its unique. It is the only free-standing pneumatic $HeliO_2$ blending device on the market today. The ability to turn the bleed control collar ON or OFF provides \pm 3% F₁₀₂ accuracy as well as quieter operation and maximal Heliox tank durations.

Part #	Flow	%	Price
5570	High Flow	70 - 30	1795.00
5580	High Flow	80 - 20	1795.00
5470	Low Flow	70 - 30	1795.00
5480	Low Flow	80 - 20	1795.00

RC MEDICAL SELLS A WIDE VARIETY OF RESPIRATORY SUPPLIES AND EQUIPMENT.
PLEASE VISIT OUR WEBSITE AT:
WWW.RCMEDICAL.COM

PRECISION MEDICAL EASY AIR COMPRESSOR



EasyAir Compressor is a high performance portable medical air compressor designed to reliably supply compressed air 24 hours a day. EasyAir is ideal for continuous use in any clinical environment. Its sound dampened case is easy to clean, allows quiet operation and is durable enough for years of heavy use.

PM15P Pressure Gauge \$719.00 PM15F Flow Gauge \$719.00

REGULATOR AND FLOWMETER ACCESSORIES









Nut and Nipple (Christmas Tree)

Part #	Description	Price
7-4A*	One piece Green	\$34.00 cs/50
7-5A*	Two piece Green	\$47.00 cs/50
7-5C*	Two piece Yellow	\$59.00 cs/50
7-5H*	Two piece Clear	\$63.00 cs/50

Plastic - For single Patient Use (Above)*		
0159	Metal Nut/Nipple	\$ 13.00 each
0159C	Metal w/ 6" chain	\$ 18.00 each

CYLINDER WRENCHES FOR SMALL AND LARGE CYLINDERS







792-D









MCW-2BC

Part # Description Price 7-92A Plastic Small Wrench, BLACK 26.00 cs/10 7-92B Plastic Small Wrench, GREEN 26.00 cs/10 7-92D Plastic Small Wrench, PINK 26.00 cs/10 66071 Plastic Small Hand Wheel 75.00 cs/10 MCW-2B Metal Small Cylinder Wrench 48.00 cs/10 MCW-2BC Metal Small Wrench with Chain 62.00 cs/10 5080 **Hudson Small Metal Wrench** 11.00 each WR-3 Metal Wrench Small & Large 18.00 each





WR-3 5080

REPLACEMENT GAUGES FOR OXYGEN REGULATORS

2" Bottom Port 1/4" NPT







2530



11/2" Back Port 1/8" NPT





2521

Low Profile 3000 PSI



MJ-C1

Part #	Description	Price
2510	15 lpm - 2" Bottom Port	\$25.00
2530	Contents (0-4000 PSI) 2" Bottom	\$25.00
2543	100 PSI 2" Bottom Port	\$25.00
2501	15 lpm - 1½" Back Port	\$25.00
2521	Contents (0-4000 PSI) 11/2" Back	\$32.00
2500	15 lpm - 11/2" Bottom Port	\$25.00
2520	Contents (0-4000 PSI) 11/2" Bottom	\$25.00
MJ-C1	Low Profile Contents (3000 psi)	\$25.00
	AIR, CO2 GAUGES	
MG-C11	AIR, CO2 (0 - 100 PSI) 2" Bottom	\$49.00
MG-C9	AIR, CO2 (0-4000 PSI) 2" Bottom	\$49.00
MG-C13	AIR, CO2 (0-4000 PSI) 11/2" Bottom	\$49.00

LITERMETER Yoke Washer

For Small Cylinder Regulators (CGA 870)



2544-2	
--------	--

Part #	Description	Price
2544-2	aluminum w/ rubber center	\$ 3.00 each
2544-2C	aluminum w/ rubber center	\$100.00 cs/50

Every therapist should have one for checking the accuracy

of regulators and flowmeters.



Part #	LPM	Price
180-5004	0-2.5	\$43.00
180-5005	0-8	43.00
180-5006	2-15	43.00

Presorted First Class U.S. Postage PAID Milford, CT Permit No. 80

ALL IN STOCK ORDERS PLACED BY 2 PM (EST) SHIPPED THE SAME DAY.

"Y" BLOCKS

Ohmeda Style

(Ohio)



2152	O ₂	193.00
2252	Air	193.00
2352	Vac	193.00

Chemetron Style





3152	O ₂	195.00
3252	Air	195.00
3352	Vac	195.00

REUSABLE NUT AND NIPPLE ADAPTER (CHRISTMAS TREE)







Main Flow Bacterial / Viral Filter

- 22mm ID x 15mm ID/22mm OD
- Bacterial filtration efficiency 99.999+%
 Viral filtration efficiency 99.99+%
 - Dead space: 42mL
 - Flow resistance: approximately
 - 1.5 cmH₂O at 60 LPM
 - Individually packaged 50 to a case

Part Number 1605 \$198.00 per case of 50

(FOR SMALL CYLINDER REGULATORS)



2544-2 \$3.00

2544-2C \$100.00/cs. 50

ADJUSTABLE HEIGHT SMALL CYLINDER CART FOR D & E CYLINDERS



#PX-1060.....\$65.00

PHONE: (888) 431-2727

WWW.RCMEDICAL.COM

FAX: (888) 337-0443